Invia[®] Motion[™]

NEGATIVE PRESSURE WOUND THERAPY SYSTEM

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Introduction	4
Overview	5
Intended Use	7
Indications for Use	7
Contraindications	8
Warnings and Safety Instructions	9
Warnings	10
Cautions	11
Physicians Orders	11
Wound Assessment	12
Dressing Technique	12
Initial Delivery	13
Invia Motion Disposables	14
Initial Set-up	15
Charge Battery for the first time	15
Disconnect from Power Source	16
Functions	17
Preparation for Use	18
Operating Instructions	21
Power Supply and Battery Management	23
Notifications and Troubleshooting	25
Overview Notifications	28
Cleaning and Sterilization	29
Disposal	30
Warranty and Maintenance	31
Accessories	32
Technical Data	33
Signs and Symbols	70
Technical Specifications	72

Federal law restricts this device to sale by or on the order of a physician.

Introduction

The Invia Motion Negative Pressure Wound Therapy System, subsequently named "Invia Motion NPWT System", is approved exclusively for the use described in these Instructions for Use. Please read the information contained within these Instructions for Use before operation of the Invia Motion NPWT System. These Instructions for Use must be kept with the device. These Instructions for Use are a general guide for the use of the product. Medical situations must be addressed by a physician.

Indications, Contraindications, Warnings, Cautions and other information to support the safe use of the Invia Motion NPWT System are included on pages 7 – 11 of these Instructions for Use.

The Invia Motion NPWT Pump is a suction pump for Negative Pressure Wound Therapy that provides treatment status through a display and acoustic signals. The Invia Motion NPWT Pump is a single patient use pump which provides continuous or intermittent operation, multiple negative pressure selection options and a lifetime of 60 days. The Invia Motion NPWT Pump is portable and can be operated using a rechargeable battery. Acoustic and optical signals are triggered for variances from the set values as well as for faults. The Invia Motion NPWT system is intended for use in all healthcare settings, including home care, by healthcare personnel and caregivers trained on the Instructions for Use. The patient or caregiver must have the visual and hearing acuity necessary to appropriately respond to notifications with the sensory and cognitive ability to understand directions from a healthcare professional pertaining to the proper use of the pump.

Overview


The Invia Motion NPWT System and the Invia Motion Carrying Case are for single patient use only.

Buttons to Operate the Invia Motion NPWT Pump

Backlight

The backlight will be activated by pressing any button. The backlight will turn off after 10 seconds if no further interaction results.

Symbols on Display

Notifications	
	Leakage
	Canister full
	Battery full
	Battery charge condition
	Battery low (symbol is blinking)
	Battery empty
	Battery charging (bars moving)
Information	
	Shows lifetime already used in days () and hours ()
	Pump is running / On-state in intermittent mode
	Off-state in intermittent mode
	Constant mode
	Intermittent mode
	Negative pressure setting in mmHg

Intended Use

The portable Medela Invia Motion Negative Pressure Wound Therapy (NPWT) system is indicated to create an environment that promotes wound healing by secondary or tertiary (delayed primary) intention by preparing the wound bed for closure, reducing edema, promoting granulation tissue formation and perfusion, and by removing exudates and infectious material. It is intended for use in hospitals, clinics, Long Term Care (LTC) and Home Care (HC) settings on adult patients with chronic, acute, subacute, traumatic, dehisced wounds, partial-thickness burns, ulcers (such as diabetic, neuropathic, pressure or venous insufficiency), flaps and grafts.

Indications for Use

The Invia Motion NPWT System is indicated for the use with the following types of wounds:

- Diabetic/Neuropathic ulcers
- Pressure or venous insufficiency ulcers
- Chronic wounds
- Traumatic wounds
- Acute or subacute wounds
- Dehisced wounds
- Partial thickness burns
- Flaps and Grafts

Intended User

The Invia Motion NPWT System is intended to be used by healthcare professionals or adequately trained lay users.

Healthcare professionals are responsible to train lay users according to the Patients Instructions for Use and explain all related safety information.

Intended Patient Population

The Invia Motion NPWT System is intended to be used on adult patients only exhibiting conditions as described in the Indications for Use.

Contraindications

The Invia Motion NPWT System is contraindicated in the presence of:

- Malignancy in the wound
- Untreated osteomyelitis
- Non-enteric and unexplored fistulas
- Necrotic tissue with eschar present
- Exposed vasculature
- Exposed nerves
- Exposed anastomotic site of blood vessels or bypasses
- Exposed organs

Warnings and Safety Instructions

WARNINGS

Indicates a potentially hazardous situation which, if not avoided, could result in death or serious injury.

CAUTIONS

Indicates a potentially hazardous situation which, if not avoided, could result in minor or moderate injury.

Safety Related Tip

Indicating useful information about the safe use of the device.

The Invia Motion Negative Pressure Wound Therapy System (subsequently named Invia Motion NPWT System) is intended for use as described in these Instructions for Use.

Medela is only responsible for the effect on BASIC SAFETY, reliability and performance of the Invia Motion NPWT System if it is used in accordance with the Instructions for Use.

Please read and observe these warnings and safety instructions before operation. These Instructions for Use must be kept with the device.

Please note that these Instructions for Use are a general guide for the use of the product. Medical situations must be addressed by a physician.

Warnings

WARNINGS

- Do not modify this equipment without authorization from the manufacturer.
- Consult the Indications For Use, Contraindications, Warnings and Cautions before using the Invia Motion NPWT Pump as a negative pressure source with the Invia Wound Dressing or Avance Dressing Kits. Failure to read and follow all instructions in this manual prior to use may result in serious or fatal injury of the patient.
- Invia Wound Therapy Instructions advise 24 hours therapy without interruption. If therapy is discontinued for more than 2 hours using foam or 24 hours using gauze, the dressing should be replaced and therapy restarted by a healthcare professional.
- Unplug the pump before cleaning.
- Do not use accessories other than those specified or sold by the manufacturer as replacement parts for internal components as it may result in increased electromagnetic emissions or decreased electromagnetic immunity of the Invia Motion NPWT Pump.
- The Invia Motion NPWT System should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the Invia Motion NPWT System should be observed to verify normal operation in the configuration in which it will be used.
- The Invia Motion NPWT Pump is intended for “Single Patient Use Only”. Reuse on another patient might cause cross contamination.
- Patients must be closely monitored for bleeding. If sudden or increased bleeding is observed, immediately stop use of the pump, apply pressure and seek immediate emergency medical attention.
- Objective indications or signs of a possible infection or complication must be addressed immediately (e.g. fever, pain, redness, increased warmth, swelling or purulent discharge). Non-observance can lead to considerable danger to the patient.
- When using the Invia Motion NPWT System in close proximity to blood vessels, organs, nerves, tendons or ligaments a protective barrier, such as a non-adherent contact layer, must be used.
- Serious or fatal injury can result from bone fragments or sharp edges (e.g. staples or hardware) that could puncture protective barriers, vessels or organs.
- Patients with spinal cord injuries (stimulation of sympathetic nervous system) who experience Autonomic Hyperreflexia during therapy should have NPWT discontinued.
- In the event that defibrillation is required, disconnect the pump from the wound dressing before the patient is defibrillated.
- Special care is advised for dressing placement and removal in order to avoid situations such as unintentional gauze or foam retention.
- This device has not been studied in pediatric patients.
- Do not bring the device into the following areas: Positron Emission Tomography (PET) or Hyperbaric Oxygen Chamber (HBO).
- The Invia Motion NPWT System is Magnetic Resonance (MR) Unsafe. Do not take the Invia Motion NPWT System into the MR environment.

Cautions

CAUTIONS

- Incorrect use may cause pain and injury to the patient. Excessive negative pressure or an infection of the wound may cause pain to the patient. In either case, the dressing must be changed and the wound inspected in detail.
- The patient should be monitored regularly according to physician instructions and facility guidelines to monitor patient comfort, therapy compliance and signs of wound infection.
- Do not use an Invia Motion Canister/Tubing Set if the sterile packaging is damaged.
- The use of mobile telephones, local area networks (LAN) including wireless (WLAN), walkie-talkies (two-way radios) and cordless telephones may affect the Invia Motion NPWT Pump. A safe distance of 3.3 feet (1 m) from the Invia Motion NPWT Pump to the device is recommended.
- The Invia Motion NPWT Pump must remain in an upright (carrying case) or lying (display side up) position during use.

The following statements describe medical conditions that may require special care for the safe and effective use of the Invia Motion NPWT System.

- Patients at high risk for bleeding and hemorrhage, including patients experiencing active bleeding or difficult wound hemostasis.
- Patients taking anticoagulants or platelet aggregation inhibitors.
- Patients with a history of vascular anastomosis or friable, irradiated, sutured or infected blood vessels.
- Patients with infected wound or osteomyelitis.
- Wounds that involve an enteric fistula.
- Use near vagus nerve (bradycardia).
- Circumferential dressing application.
- The Negative Pressure Wound Therapy must be used 24 hours per day without interruption. If the pump is stopped for more than the time frames shown below, the dressing must be changed and therapy restarted.

Gauze dressing: Change dressing if the pump has stopped for more than 24 hours.

Foam dressing: Change dressing if the pump has stopped for more than 2 hours.

- Consider the patient's size and weight when prescribing this device.
- Consider mode of therapy – intermittent versus continuous.

Physicians Orders

A medical professional must assess each wound when applying this therapy and determining negative pressure intensity according to the wound characteristics. Therapy or setting changes may only be done under the order of a physician.

Wound Assessment

CAUTION

Patient Monitoring: The patient should be monitored regularly according to the physician's instructions and facility guidelines to check for patient comfort, therapy compliance and signs of infection.

WARNING

Objective indications or signs of a possible infection or complication must be addressed immediately (e.g. fever, pain, redness, increased warmth, swelling or purulent discharge). Non-observance can lead to considerable danger to the patient.

Observe the wound/periwound tissue and exudate for signs of infection or other complications. The most common signs of infection include redness, tenderness, fever, swelling, itching, increased warmth in the wound area, strong odor or purulent discharge. Additional symptoms include nausea, vomiting, diarrhea, headache, dizziness, fainting, sore throat with swelling of the mucous membranes, disorientation, high fever ($>102^{\circ}\text{F}$, 38.8°C), refractory hypotension, orthostatic hypotension, or erythroedema (a sunburn-like rash). More serious complications of infection include pain, discomfort, fever, gangrene, toxic or septic shock. If more serious complications of infection occur, discontinue therapy and consult a healthcare professional immediately.

Dressing Technique

Consult the appropriate Invia Wound Dressing or Avance Dressing Kit Instructions for Use for information regarding dressing applications. Perform a thorough wound cleansing per physician orders prior to dressing applications.

Routine dressing changes should occur every 48–72 hours. Dressing changes for infected wounds should be considered more frequently.

Initial Delivery

Upon delivery, check the package containing the Invia Motion NPWT System for completeness and general condition.

REF 087.4000 Invia Motion Negative Pressure Wound Therapy System

- | | |
|---------------------|--|
| REF 200.4205 | 1 Invia Motion NPWT Pump |
| REF 087.4001 | 1 Invia Motion Canister/Tubing Set 0.15l (sterile, single use) |
| REF 200.4212 | 1 Invia Motion Power Supply |
| REF 087.4003 | 1 Invia Motion Carrying Case |
| REF 200.4168 | 1 Invia Motion Negative Pressure Wound Therapy System Clinician Instructions for Use |
| REF 200.4169 | 1 Invia Motion Negative Pressure Wound Therapy System Patient Instructions for Use |

☞ When first turned on, the Invia Motion pump needs to be attached to a wall outlet in order to fully charge the battery.

Invia Motion Disposables

Invia Motion Canister/Tubing Set 0.15l (150 ml), with Solidifier (Sterile, Single Use)

Capacity: 0.15l (150 ml)

Material: Polypropylene, Silicon

Sterilization Method: Ethylene Oxide


Intended Use

The Medela Invia Motion Canister/Tubing Set 0.15l (150 ml) is an accessory of the Medela Invia Motion NPWT System and it is intended to collect removed wound exudate and infectious material.


Initial Set-up

Attach the US Plug to the Invia Motion Power Supply

Charge Battery for the first time

Plug the Invia Motion Power Supply into a wall outlet and open the small protective cover of the pump for the Direct Current (DC) socket and attach the power supply. Moving bars on the battery display will indicate charging is occurring. [].

 Check the charging state of the battery periodically.

1. Battery empty


2. Battery is being charged
(moving bars)

3. Battery is fully charged

Disconnect from Power Source

Disconnecting the Power Source

When the battery is fully charged (after approximately 4 hours) the following non-blinking symbol  will appear on the display. To disconnect the Invia Motion NPWT Pump from the electrical source, hold the end of the plug and remove it from the pump.

Close the DC socket by pressing the small protective cover over the opening.

Functions

Pressure Regulation

The pressure is checked periodically by the Invia Motion NPWT Pump and maintained according to the settings entered by a healthcare professional.

Canister Safety Chamber

CAUTION

The Invia Motion NPWT Pump must remain in an upright (carrying case) or lying (display side up) position during use.

The Invia Motion NPWT Pump is protected against the penetration from solid/fluid substances by a hydrophilic bacteria filter. If this filter fails, the Invia Motion NPWT Pump must be replaced.

In the event of inappropriate positioning of the pump, the special construction of the safety chamber in the upper region of the canister protects the hydrophilic filter from premature clogging to maintain the pressure.

If the canister is full and the Invia Motion NPWT Pump is placed in an inappropriate position, fluids will flow into the safety chamber and clog up the hydrophilic overflow protection/bacteria filter. In this situation, the pressure is interrupted and an acoustic information signal will sound. The Invia Motion Canister/Tubing Set must be replaced.


Preparation for Use

 To be used by healthcare professionals or adequately trained caregivers only.

Checklist

Prior to using the Invia Motion NPWT System, it is important to check for the following:

- Completeness and general condition of the packaging
- Damage to the power cord and plug
- Damage to the pump
- Proper functioning of the device
- Battery status

1. Open sterile packaging.

CAUTION


Do not use the Invia Motion Canister/Tubing Set if the sterile packaging is damaged.

2. Insert the Invia Motion Canister/Tubing Set into the Invia Motion NPWT Pump connecting port.

3. Rotate the Canister/Tubing Set to the right to secure.

4. Connect the Invia Motion Tubing Connector to the wound drain tubing or the tubing of the External Suction Interface (ESI).


For dressing applications, please refer to the Instructions for Use provided with the Invia and Avance Dressing Kits.


Switch Invia Motion NPWT Pump On

1. Press [] to switch the Invia Motion NPWT Pump on.

2. When turned on for the first time, the Invia Motion NPWT Pump will start up in the default factory setting (constant mode with pressure level of -125mmHg).


If the Invia Motion NPWT Pump is turned off, the setting will be as previously set when the pump is turned back on.

Lifetime, Countdown and Safety Time

The Invia Motion NPWT Pump has a lifetime of 60 days. After one hour of use, the countdown of 60 days is activated and will continue even when the pump is turned off. During the first hour of use, the Invia Motion Pump can be turned on and off without starting the countdown.

After 60 days, a blinking symbol on the display [] will indicate the expiration of the Invia Motion NPWT Pump.

- ☞ The Invia Motion NPWT Pump provides an additional safety usage time of 72 hours. During this time the device can be operated as usual. An acoustic signal will occur after 72 hours and the negative pressure will be maintained for one hour. The acoustic signal cannot be muted and the pump will shut down automatically after one hour unless it is turned off earlier. Once turned off, the Invia Motion Pump cannot be turned on again and has to be disposed of according to local regulations.

Setup Carrying Case

The Invia Motion NPWT Pump is intended to be worn in a carrying case during operation. An Invia Motion Carrying Case is provided with the Invia Motion NPWT Pump.

Attach shoulder/belt strap onto the carrying case as a first step.

1. Open the carrying case.
2. Slide the Invia Motion NPWT Pump into the designated upper pocket and close the Velcro strap.
3. Place the excess tubing into the designated lower pocket and close the Velcro strap.
4. Fold the pump and tubing sections of the carrying case together and align the inner Velcro tabs to secure.
5. Secure the bag by closing the lid and fastening the button.

The carrying case strap can be worn over the shoulder or around the waist.

To wear over the shoulder, attach the strap clips to the carrying case. To wear as a belt, thread the strap through the two loops located on the back of the case.

Shoulder strap

Belt strap

Operating Instructions

Invia Motion NPWT Pump Settings

Change Pressure Level

Press [] and hold for 3 seconds to change the pressure level. The set pressure value symbol on the display will start to blink. You now have 5 seconds to press [] for a new pressure setting, otherwise the pump remains at the previous setting. Each time the button is pressed the next available pressure will appear. When the appropriate pressure level appears on display, the pump will adopt it automatically after 5 seconds. The pressure may only be changed if instructed by a physician or qualified healthcare professional.

Available pressure levels

mmHg	-60	-80	-100	-125	-150	-175

Change Therapy Mode


The Invia Motion NPWT Pump is able to deliver constant or intermittent (5 minutes on/2 minutes off) suction.


Constant mode

Intermittent mode

Press [] for 3 seconds to change the therapy mode. The therapy mode symbol on the display will start to blink. You now have 5 seconds to press [] to change the mode. Otherwise the pump remains at the previous mode. Each time the button is pressed the next available therapy mode will appear. When the appropriate mode symbol appears, the pump will adopt it automatically after 5 seconds. The therapy mode may only be changed if instructed by a physician or a qualified healthcare professional.

Pump in constant mode

Pump in intermittent mode

Mute / Information

Press [] to obtain information about the operation time used or to mute the acoustic signal.


Change Invia Motion Canister/Tubing Set

The Invia Motion Canister/Tubing Set has to be changed on the basis of a visual check or according to the instructions on the display.

1. Wash hands.
2. Switch the Invia Motion NPWT Pump off by pressing [] for 3 seconds.
Press only once.
3. Disconnect the Invia Motion Tubing Connector from the wound or External Suction Interface (ESI) drain tubing connector. Have a paper towel or tissue available to cover the tip of the used tubing to protect against wound fluid leakage.
4. Remove the pump from the carrying case, detach the full Canister/Tubing Set by rotating the canister to the left. Discard the used Canister/Tubing Set appropriately.

5. Open the sterile packaging of a new Invia Motion Canister/Tubing Set.
6. Remove the Canister/Tubing Set from the packaging and attach it to Invia Motion NPWT Pump.

For dressing applications,
please refer to Instructions
for Use provided with Invia
and Avance Dressing Kits.

7. Connect Invia Motion Tubing Connector to the wound drain tubing.
8. Switch the Invia Motion NPWT Pump on by pressing [] and slide it back into the carrying case.

Power Supply and Battery Management

The Invia Motion NPWT Pump is operated by an internal rechargeable NiMH battery and while attached to a power source. While in use and connected to a power source, the battery is recharged. The available battery time from a single charge depends upon the run-time of the pump. This is influenced by the size of the wound, air leaks in the dressing, and the set pressure. The pump only runs if the pressure drops below the set pressure. If the Invia Motion NPWT Pump is running continuously, a maximum of 4 hours of battery time is expected before it is necessary to be recharged. The expected battery time generally exceeds 10 hours. If the Invia Motion NPWT Pump is not in use, the battery must be charged approximately once every 6 months to ensure optimal functioning.

 Before attaching the Invia Motion NPWT Pump to a power source, check to ensure your local power supply is the same as the voltage indicated on the specification plate located on the back of the pump.

Battery Low

Approximately 30 minutes before the battery is fully empty, an acoustic signal will sound and the battery symbol [] will start to blink. The acoustic signal can be muted, but the battery symbol continues to blink until the battery is empty. The set pressure is maintained during the 30 minutes, but the battery should be recharged as soon as possible.

Battery Empty

The battery symbol will stop blinking and an acoustic signal will sound for 10 minutes prior to the battery fully discharging, unless the pump is switched off earlier by pressing [] and holding for 3 seconds. The acoustic sound cannot be muted during the 10 minutes. If the Invia Motion NPWT Pump is not connected to a power source the pump will switch off automatically after 10 minutes. The negative pressure is not maintained after the pump turns off.

Connecting Power Source

Plug the Invia Motion Power supply into a wall outlet and open the small protective cover to access the DC socket. Plug the coaxial power supply into the DC socket. Moving bars on the battery display will indicate charging is occurring [].

1. Battery empty

2. Battery is being charged
(moving bars)

3. Battery is fully charged

Disconnecting Power Source

When the battery is fully charged (after approximately 4 hours) the following non-blinking symbol  will appear on the display. To disconnect the Invia Motion NPWT Pump from the electrical source, hold the end of the plug and remove it from the pump.

Close the DC socket by pressing the small protective cover over the opening.

Notifications and Troubleshooting

WARNING

Negative Pressure Wound Therapy should run 24 hours a day without interruption. If therapy is discontinued for more than: 2 hours with a foam dressing or 24 hours with a gauze dressing, the dressing should be replaced and therapy restarted by a healthcare professional.

If the Invia Motion NPWT Pump detects an irregularity in the operation of the pump, an acoustic signal sounds and a notification will appear on the display.

1. An acoustic signal sounds.
2. Press [] to mute the acoustic signal for 60 seconds (the button can be pressed every time the signal sounds, with the exception of "Battery Empty" and after the countdown of 60 days plus 72 hours of safety time is used up).


Every acoustic signal can be permanently muted by switching Invia Motion NPWT Pump off (press [] and hold for 3 seconds). Turn the pump back on to continue treatment.

3. If the Invia Motion NPWT Pump is in the carrying case, open the lid of the carrying case and release the Velcro.

4. Check the display for one of the following visual symbols.

Canister Full


Troubleshooting:

Change the Invia Motion Canister/Tubing Set.


Leak in System

Troubleshooting:

Check the dressing and Canister/Tubing Set connections for leaks. If a leak is identified at the dressing site, press down on dressing to seal.

Battery Low


Troubleshooting:

Charge the battery.

Approximately 30 minutes before the battery is fully empty, an acoustic signal will sound and the battery symbol [] will start to blink. The acoustic signal can be muted, but the battery symbol continues to blink until the battery is empty. The set pressure is maintained during the 30 minutes, but the battery should be recharged as soon as possible.


Battery Empty

Troubleshooting:

Charge battery

The battery symbol will stop blinking and an acoustic signal will sound for 10 minutes prior to the battery fully discharging, unless the pump is switched off earlier by pressing [] and holding for 3 seconds. The acoustic signal cannot be muted during the 10 minutes. If the Invia Motion NPWT Pump is not connected to a power source the pump will switch off automatically after 10 minutes. The negative pressure is not maintained after the pump turns off.

Info Used Pump Time

The Invia Motion NPWT Pump has a limited lifetime of 60 days. After one hour of use, the countdown of 60 days is activated. During the first hour of use, the Invia Motion NPWT Pump can be turned on and off without starting the countdown. The used time within the first hour is summed up and the countdown is triggered at 60 minutes of usage. Press [] for information about the operation time used.

The example below displays 6 days and 23 hours of operation time used. The remaining lifetime of Invia Motion NPWT Pump would be 53 days and 1 hour.

Countdown Expired

After 60 days, a blinking symbol on the display [] will indicate the expiration of the Invia Motion NPWT Pump.

The Invia Motion NPWT Pump provides an additional safety usage time of 72 hours. During this time the device can be operated as usual. An acoustic signal will sound after 72 hours and the negative pressure will be maintained for one hour. The acoustic signal cannot be muted and the pump will shut down automatically after one hour unless it is turned off earlier. Once turned off, the Invia Motion NPWT Pump cannot be turned on again and has to be disposed of according to local regulations.

Overview Notifications

	Visual Symbol on Display	Notification Description	Acoustic Signal	Troubleshooting	Remarks/Potential Cause
Notification	 + acoustic signal!	Canister full	Can be muted for 60 seconds by pressing. 	Change the Invia Motion Canister/ Tubing Set.	The filter can get wet by improper positioning of the pump (upside down or display side down).
	 + acoustic signal!	Leak in system	Can be muted for 60 seconds by pressing. 	Check the system for leaks.	Dressing or Canister/ Tubing Set connections are not properly sealed.
	 Blinking! + acoustic signal!	Battery low	Can be muted by pressing. 	Charge battery.	The visual symbol blinks all the time until battery is empty. The remaining time of the battery is approximately 30 minutes.
	 + acoustic signal!	Battery empty	Can be muted by switching off the pump.	Charge battery.	The remaining time of the battery is approximately 10 minutes.
	 Blinking!	Countdown expired	The blinking signal can be suppressed by pressing. 	Dispose of the Invia Motion Pump. Attach another device if necessary.	The Invia Motion NPWT Pump provides an additional safety usage time of 72 hours. During this time the pump can be operated as usual. A signal will occur after 72 hours.
Internal fault	 + acoustic signal!	Internal fault	Can be muted by switching off the pump.	Switch the Invia Motion NPWT Pump off, release and refix canister. Turn the pump back on. Contact Medela customer service.	If the fault message repeats, inform Medela customer service.

Contact your local Medela customer service representative for assistance with product operations.

Cleaning and Sterilization

Cleaning

WARNING

Unplug the pump before cleaning.

The Invia Motion NPWT Pump and carrying case are for single patient use only. Occasional cleaning may be necessary.

Wipe both with a damp cloth or use a non-abrasive detergent.

Sterilization

The Invia Motion NPWT Pump and Accessories cannot be sterilized. The Canister/Tubing Set is a disposable item and may not be reused or resterilized. Replace after every use!

The Invia Motion Canister/Tubing Set is a single patient use product and is not intended to be reprocessed. Reprocessing may cause the loss of mechanical, chemical and/or biological characteristics.

WARNING

Reuse on another patient might cause cross contamination.

General notes

Comply with facility or institutional guidelines on hygiene and infection control.

Dispose of fluids such as blood and fluids, along with all parts contaminated with them in accordance with facility requirements. Wear suitable gloves for cleaning/disinfection and disposal.

Disposal

Turn the Invia Motion NPWT Pump off by pressing [] and hold for 3 seconds.

Disconnect the Canister/Tubing Set from the pump.

Dispose of the Canister/Tubing Set and pump separately and according to local regulations.

Take care to dispose of the Invia Motion NPWT System and its accessories in accordance with your local disposal regulations.

Disposal of Electronic Equipment

The symbol above means that the electrical and electronic equipment must not be disposed as normal household refuse. The correct disposal of this device protects and prevents possible damage to the environment or human health. For more information about proper disposal, contact Medela or your local healthcare professional. This symbol is only valid in the European Union. Please adhere to relevant state laws and regulations for the disposal of electrical and electronic equipment.

Warranty and Maintenance

Warranty

The Invia Motion NPWT System warranty period is for 2 years after the date of delivery, but expires after 60 days of usage plus 72 hours of safety time, if used in accordance with the Instructions for Use. The manufacturer is not liable for any damage or consequential damage caused by incorrect operation, inappropriate usage or unauthorized persons using the device.

Maintenance

The Invia Motion NPWT Pump is maintenance free and will not require service.

If an Invia Motion NPWT Pump fails within the warranty period due to a manufacturing defect, the pump will be replaced. The original pump will need to be returned to the supplier.

Safety-Related Checks

The Invia Motion NPWT System is a device in electrical protection Class II. The safety-related checks are confined to visual inspection of the pump housing and power supply for damage. This check must be performed prior to each use.

Devices of protection Class II do not have a protective earth conductor. There is no need to check the earth leakage current.

The Invia Motion NPWT Pump enclosure is made entirely of insulated material. Tests of the enclosure leakage current using common measuring instruments will therefore not reveal measurable values.

The Invia Motion NPWT Pump does not have patient circuits or functional earth connections.

Accessories

The Invia Motion NPWT System was verified in combination with the accessories listed below. For correct and safe operation use the Invia Motion NPWT System with these accessories only:

087.4000 Invia Motion Negative Pressure Wound Therapy System

- 1 Invia Motion NPWT Pump
- 1 Invia Motion Canister/Tubing Set 0.15l (sterile, single use)
- 1 Invia Motion Carrying Case
- 1 Invia Motion Power Supply
- 1 Invia Motion Negative Pressure Wound Therapy System Clinician Instructions for Use
- 1 Invia Motion Negative Pressure Wound Therapy System Patient Instructions for Use

087.4001 Invia Motion Canister/Tubing Set 0.15l

- 10 Invia Motion Canister/Tubing Set 0.15l (sterile, single use)

Wound Dressing Kits

Avance Foam

- REF** 662151 Avance Foam Dressing Kit incl. Transparent film, Transfer Pad - Small
- REF** 662251 Avance Foam Dressing Kit incl. Transparent film, Transfer Pad - Medium
- REF** 662351 Avance Foam Dressing Kit incl. Transparent film, Transfer Pad - Large

Invia Gauze

- REF** 877061 Tier 1 Dressing Set 10mm Flat Drain
- REF** 877063 Tier 1 Dressing Set 15Fr Round Channel Drain
- REF** 877070 Tier 1 Dressing Set ESI
- REF** 877069 Tier 1 Large Dressing Set 10mm Flat Drain
- REF** 877074 Tier 1 Large Dressing Set 15Fr Round Channel Drain
- REF** 877064 Tier 1 Large Dressing Set 19Fr Round Perforated Drain
- REF** 877071 Tier 1 Large Dressing Set ESI

The Invia Gauze Dressing kits are only available in the US. The usage of Avance Foam Dressing kits in combination with the Invia Liberty or Invia Motion is only approved for the US. Invia Gauze Dressing kits and Avance Foam Dressing kits provide their own Instruction for Use. All Invia Gauze Dressing kits and Avance Foam Dressing kits have a separate registration approval for the US market.

EMC

The Invia Motion NPWT System is EMC-tested in conformity with the requirements of IEC 60601-1-2:2007 and can be used in the vicinity of other EMC-tested devices that fulfill the requirements of the relevant IEC 60601-1-2 standard. Untested HF (high-frequency) sources, radio networks or the like can influence the operation of the device and may not be operated in combination with the system.

Essential Performance: Medela has determined that the Invia Motion NPWT System does not have any functions, the failure of which to work properly, could impact the safety of the system. As such, there is no essential performance as defined in the IEC 60601-1-2 standard. However, the Invia Motion was tested per IEC 60601-1-2 ensuring that the pump could maintain the correct vacuum during testing. Disturbances higher than those specified in IEC 60601-1-2 may cause insufficient battery charging or may lead to an internal fault. See the chapter Overview Notifications for actions to take in the case of an internal fault.

Certain mobile phones may interfere with the Invia Motion NPWT system when receiving and making calls close to the pump. Prior to using a mobile phone near the Invia Motion NPWT system, check that it continues to function while making or receiving a call. If you observe a change in the vacuum level or an error condition, move the Invia Motion NPWT system at least 1 m (3.3 ft) from the mobile phone when making or receiving calls. Turn the Invia Motion NPWT system on and off again to reset it after an error condition. Invia Wound Therapy Instructions advise 24 hours of therapy without interruption. If therapy is discontinued for more than 2 hours using foam or

24 hours using gauze, the dressing should be replaced and therapy restarted by a healthcare professional.

The Invia Motion NPWT Pump is a medical device that requires special safety precautions and must be installed and placed in operation in accordance with the attached EMC information. Portable and mobile RF communication devices (mobile telephones) can affect Invia Motion NPWT Pump.

WARNING

Do not use other accessories than those specified or sold by the manufacturer as replacement parts for internal components as it may result in increased emissions or decreased immunity of the Invia Motion NPWT Pump.

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IEC 60601-1-2, Table 1

Electromagnetic Emissions		
The Invia Motion NPWT System is intended for use in the electromagnetic environment specified below. The customer or the user of Invia Motion NPWT System should assure that they are used in such environment.		
Emission Tests	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The Invia Motion NPWT System uses RF energy only for their internal function. Therefore, its RF emissions are very low and are not likely to cause any interference with nearby electronic equipment.
RF emissions CISPR 11	Class B	
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations / flicker emissions IEC 61000-3-3	Complies	The Invia Motion NPWT System is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.

WARNING

The Invia Motion NPWT System should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the Invia Motion NPWT System should be observed to verify normal operation in the configuration in which it will be used.

IEC 60601-1-2, Table 2

Electromagnetic Immunity

The Invia Motion NPWT System is intended for use in the electromagnetic environment specified below. The customer or the user of Invia Motion NPWT System should assure that they are used in such environment.

Immunity Tests	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic Discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 16.5 kV air	The relative humidity should be at least 5 %.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input / output lines	± 2 kV for power supply lines not applicable	Mains power quality should be that of a typical domestic or hospital environment.
Surge IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV differential mode not applicable	Mains power quality should be that of a typical domestic or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	$<5\% U_T$ ($>95\%$ dip in U_T) for 0.5 cycle $40\% U_T$ (60% dip in U_T) for 5 cycles $70\% U_T$ (30% dip in U_T) for 25 cycles $<5\% U_T$ ($>95\%$ dip in U_T) for 5 sec	$<5\% U_T$ ($>95\%$ dip in U_T) for 0.5 cycle $40\% U_T$ (60% dip in U_T) for 5 cycles $70\% U_T$ (30% dip in U_T) for 25 cycles $<5\% U_T$ ($>95\%$ dip in U_T) for 5 sec	Mains power quality should be that of a typical domestic or hospital environment. If the user of the Invia Motion NPWT System requires continued operation during power mains interruptions, it is recommended the Invia Motion NPWT System be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	30 A/m	Separation distance $r = \frac{I}{188}$ Where I is the current in amperes in a power bus or an appliance wire and r is the separation distance in meters.

NOTE: U_T is the a.c. mains voltage prior to application of the test level.

Electromagnetic Immunity The Invia Motion NPWT System is intended for use in the electromagnetic environment specified below. Users of the Invia Motion NPWT System should assure that they are used in such environment.			
Immunity Tests	IEC 60601 test level	Compliance level	Electromagnetic environment-guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms	Portable and mobile RF communications equipment should be used no closer to any part of the Invia Motion NPWT System, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = 1.2 \sqrt{P}$ $d = 0.2 \sqrt{P}$ 80 MHz to 800 MHz $d = 1.2 \sqrt{P}$ 800 MHz to 2.5 GHz Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m) Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,^a should be less than the compliance level in each frequency range.^b Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF IEC 61000-3	3 V/m 80 MHz to 2.5 GHz	20 V/m 80 MHz to 800 MHz 6 V/m 800 MHz to 2.5 GHz	
NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies. NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people. ^a Field strengths from fixed RF transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which Invia Motion NPWT System are used exceeds the applicable RF compliance level above, Invia Motion NPWT System should be observed to verify normal operation. If abnormal operation is observed, additional measures may be necessary, such as reorienting or relocating Invia Motion NPWT System. ^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.			

Recommended separation distance between portable and mobile RF communications equipment and The Invia Motion NPWT System

The Invia Motion NPWT System is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. Users of the Invia Motion NPWT System can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Invia Motion NPWT System as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter	Separation distance according to frequency of transmitter		
	m		
W	150 kHz to 80 MHz $d = 1.2 \sqrt{P}$	80 MHz to 800 MHz $d = 0.2 \sqrt{P}$	800 MHz to 2.5 GHz $d = 1.2 \sqrt{P}$
0.01	0.12	0.02	0.12
0.1	0.38	0.07	0.38
1	1.2	0.2	1.2
10	3.8	0.7	3.8
100	12	2	12

For transmitters rated at a maximum output power not listed above, the recommended separation distance (d) in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Signs and Symbols / Signos y símbolos

This Symbol indicates to consult Instructions for Use.

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This Symbol indicates Safety related Tip.

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This Symbol indicates a CAUTION or WARNING associated with the Device (see page 3).

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This Symbol indicates to not dispose the Device together with Household Refuse (for EU only).

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This Symbol indicates the Date of Manufacture (four Digits for the Year and two Digits for the Month).

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This Symbol indicates the Name and the Address of the Manufacturer.

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This Symbol indicates the Device is sterilized using ethylene oxide.

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This Symbol indicates a Prescription Device. Federal Law restricts this Device to sale by or on the order of a licensed healthcare practitioner. (for US only).

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This Symbol indicates a type BF applied Part.

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This Symbol indicates manufacturer's Catalog Number.

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This Symbol indicates manufacturer's Serial Number.

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This Symbol indicates manufacturer's Batch Code.

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This Symbol indicates the Protection against the ingress of dripping water and solid foreign objects of 12,5 mm and greater.

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This Symbol indicates that the Device should not be used after the End of the Year and Month shown.

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Medela and Invia are registered trademarks of Medela Holding AG.
Avance is a registered trademark of Mölnlycke.

CE 0123

If you need assistance contact:

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Medela AG
Medela Inc. USA