Use and Care Manual









AUTO
ICON^{AT}








PREMO
ICON^{PT}










NOVO
ICON^{NT}

A

B

C

English	A
简体中文版 (Chinese Simplified)	B
Dansk	C
Deutsch	D
Español	E
Français	F
Greek	G
Italiano	H
日本語 (Japanese)	I
한국어 (Korean)	J
Nederlands	K
Norsk	L
Português	M
Suomi	N
Svenska	O
繁體中文版 (Chinese Traditional)	P
Turkish	Q

TABLE OF CONTENTS

1.	Before you Start	A - 1
2.	Symbol Definitions	A - 1
3.	General Warnings and Cautions	A - 2
4.	Intended Use	A - 3
5.	Quick Start Guide	A - 4
6.	Getting to Know Your F&P ICON™	A - 5
	Replacement Parts	A - 6
7.	Menu System	A - 7
	(1) On/Off and Ramp	A - 7
	(2) Humidity	A - 7
	(3) Alarm On/Off	A - 7
	(4) Alarm Time	A - 7
	(5) Alarm Volume	A - 7
	(6) Clock Time	A - 7
	(7) Brightness	A - 7
	(8) Sleep Data	A - 8
	(9) SensAwake™ On/Off or Altitude Adjustment	A - 8
	(10) ECO Mode	A - 8
	(11) Display Pressure	A - 8
	(12) Service Setting	A - 8
8.	Communicating Sleep Data	A - 9
	Telephone Reporting	A - 9
	Downloading Data to the SmartStick™	A - 9
9.	Cleaning and Maintenance	A - 9
10.	Product Specifications	A - 10
11.	Device and Consumables Disposal Instructions.....	A - 11
12.	Troubleshooting	A - 11
13.	F&P ICON™ Models and Features Matrix	A - 11

1. BEFORE YOU START

Please read this manual carefully before use and save these instructions for future reference.

If you feel that your device or any accessories are not operating correctly, please contact your healthcare provider.

NOTE: This device contains no serviceable parts. Refer all therapy and equipment queries to your healthcare provider.

IMPORTANT CONTACT DETAILS

Home Medical Equipment Provider

Doctor/Physician

Health Insurance Provider

CAUTION: USA Federal law restricts this device for sale by or on the order of a physician.

2. SYMBOL DEFINITIONS

Type BF Applied Part

Alternating Current

ATTENTION
Consult accompanying documents

Class II Double-insulated

Do not throw away

CE 0123

Conforms with medical device directive 93/42/EEC

IPX2

Drip-proof

Date of manufacture

REF

Catalogue number

SN

Serial number

Rx Only

Prescription only

EC REP

Authorized representative in the European community

Z078 NZ Radio Interference C-tick mark

3. GENERAL WARNINGS AND CAUTIONS

This manual refers to the **F&P ICON™** series model as “the device”.

The device treats Obstructive Sleep Apnea (OSA) by delivering a flow of continuous positive airway pressure (CPAP) at a level prescribed by the physician, to splint open the airway and prevent airway collapse.

CONTRAINDICATIONS

- Research indicates that the following pre-existing conditions may contraindicate the use of positive pressure for some patients: pneumothorax, bullous lung disease, pneumocephalus, cerebrospinal fluid leak, recent cranial surgery or trauma, abnormalities of the cribriform plate, pathologically low blood pressure or in patients whose upper airways are bypassed.

PRECAUTIONS

- The safety and effectiveness of positive pressure has not been established in patients with respiratory failure or chronic obstructive pulmonary disease.
- The safety and effectiveness of the auto-adjusting device has not been established in patients with congestive heart failure, obesity hypoventilation syndrome or central sleep apnea.

ADVERSE EFFECTS

- Nosebleeds, ear and sinus discomfort may occur from the use of positive pressure therapy.

WARNINGS

To avoid electric shock from the device:

- Only operate if the device, power cord and plug are dry and in good working order.
- If water damage occurs to your device, disconnect the power cord, discontinue use immediately and seek advice from your healthcare provider.
- Do not store or use the device where it can be pulled into water.

To avoid choking, or inhalation of a foreign object:

- Never place any non-approved objects into any opening of the device, breathing tube or mask.
- Ensure the recommended filter is fitted to the device before use.
- Ensure the breathing tube is positioned so it can not become entangled with the body or furniture during sleep.

To avoid burns:

- Do not fill the chamber with hot water.
- After use, wait for the water to cool before touching, carrying or emptying the Water Chamber.
- Do not touch the water in the chamber while the device is operating.
- Position the ThermoSmart™ Breathing Tube so it is uncovered and free from bedding or other materials.

To ensure optimal therapy:

- Do not operate the device, chamber or breathing tube if it is dropped, damaged or not working as intended.
- Pressure adjustments should only be made by a qualified healthcare provider.
- Only use masks, breathing tubes and accessories compliant with ISO 17510-2, distributed for use with this device, and recommended by Fisher & Paykel Healthcare or your healthcare provider.
- Do not block the exhaust flow holes of the mask as they are designed to allow a continuous flow of air out of the mask and CO₂ re-breathing may occur.
- In the event of power failure, machine malfunction or if the device is turned off, remove the mask immediately as the flow through the mask may be insufficient to clear all exhaled gas and accordingly CO₂ re-breathing may occur which can be hazardous.
- On models without automatic altitude adjustment, ensure the altitude level is manually adjusted to ensure optimal pressure delivery.
- Do not use the device without a filled water chamber in place.

Other:

- Refer all repairs and maintenance to your healthcare provider.
- Ensure the device is stored on a level surface, lower than head height, to prevent water entering the tubing and the device enclosure.
- The device is not intended to be used as a life-supporting device.
- The device complies with the electromagnetic compatibility requirements of IEC 60601-1-2. In certain circumstances the device may affect or be affected by nearby portable mobile radio frequency communication equipment, due to the effects of electromagnetic interference. If this should happen, try moving your device or the location of the equipment causing interference, or alternatively consult your healthcare provider.

Using supplemental oxygen with your device:

1. Supplemental oxygen can be administered at the mask end of the breathing tube.
NOTE: At a fixed flow rate of supplemental oxygen the inhaled oxygen concentration will vary, depending on the pressure settings, the patient's breathing pattern, mask selection and leak rate.
2. Ensure there is no obstruction downstream of the Oxygen Port as this can affect the delivered oxygen concentration.
3. Oxygen concentration should be measured at the point of delivery to the patient.
4. Avoid the risk of fire:
 - Only use oxygen when the device is operating. If the device is turned off it can lead to accumulation of oxygen within the device.
 - Ensure adequate ventilation is provided around the device.
 - Remove any source of ignition, such as: cigarettes, an open flame, or materials which burn or ignite easily at high oxygen concentration.
 - Keep oxygen regulators, cylinder valves, tubing, connections and all other oxygen equipment away from oil, grease or greasy substances. Spontaneous and violent ignition may occur if these substances come into contact with oxygen under pressure.

CAUTIONS**To prevent water damage to the device:**

- Remove the Water Chamber from the device before filling.
- Do not fill the Water Chamber above the maximum level.
- Do not move, carry, transport or store the device with water in the chamber.

General:

- Only use the device within the Operating Conditions specified in Section 10.
- Position the device so the power cord connection to the power supply is easily accessible.
- Only clean the device in accordance with the cleaning instructions set out in Section 9 and only when it is disconnected from the power supply.
- Only use the F&P ICON™ SmartStick™ with the device. Use of any other USB drives may cause data corruption. Do not attempt to change the directories or view the data without software distributed or designed for use with the F&P ICON™.

4. INTENDED USE

The device is for use on adult patients for the treatment of OSA. The device is for use in the home or sleep laboratory.

5. QUICK START GUIDE

Refer to **Part A** of the Diagram Sheet at the beginning of this manual for illustrative examples:

- 1 Select the appropriate plug from the Carry-Case and attach to the end of the power cord. Place the device on a firm flat surface and plug the power cord into the wall socket.
- 2 Reach inside the Water Chamber with your fingers to take hold of the circular chamber grip. Pull the Water Chamber up and out of the device.
CAUTION: Do not fill the Water Chamber while it is in the device; the Water Chamber should always be removed for filling.
- 3 Fill the Water Chamber with distilled water up to the maximum water line, which is marked on the side of the Water Chamber.
CAUTION: Do not use the device without a filled Water Chamber.
- 4 To replace the Water Chamber, the arrow on the circular chamber grip should face toward the Elbow.
- 5 Fit the Chamber Lid by lining up the large dot on the lid ● with the large dot on the top of the device ●. When the two symbols are aligned, apply downward force and turn clockwise to align the large dot with the small dot on the top of the device.
- 6 When the Chamber Lid is locked into place, it should be flush, with no gaps for air to leak. It can now be used as a handle for lifting and carrying.
NOTE: The Chamber Lid must be correctly locked for optimal pressure delivery.
- 7 The Chamber Lid can be unlocked by twisting in an anti-clockwise direction, then remove by lifting away from the device.
- 8 Attach the breathing tube to the Elbow and connect to your mask.
- 9 Press the SmartDial™ to start pressure delivery. The Home Screen will light up and the clock-face of the Menu System will rotate to indicate therapy is starting. If Ramp is required, press the SmartDial™ again and hold for 3 seconds to activate the Ramp. The Ramp symbol  will light up on the display.
- 10 Turn the SmartDial™ to select and adjust settings. Refer to Section 7 for a description of the Menu System.

Your device is now ready for use.

6. GETTING TO KNOW YOUR F&P ICON™

Description of Parts and Instructions for Use

Refer to **Part B** at the beginning of this manual for diagrams of each part. Also refer to **Section 5 Quick Start Guide** for additional information.

1 Water Chamber

- Remove the Water Chamber from the device before filling. Fill to the maximum water line (420 ml), using distilled water only. Do not use hot water.

CAUTION: Do not fill the Water Chamber while it is inside the device and do not fill above the maximum water line.

- To replace the Water Chamber into the device, the arrow on the circular chamber grip should face toward the elbow.

NOTE: Replace water before each use.

2 SmartDial™

The SmartDial™ is for operating the device and navigating the Menu System. The SmartDial™ works on a simple *Press, Turn, Press* principle. For example:

- Press* the SmartDial™ to turn the device on or off.
- Turn* the SmartDial™ to find the required setting.
- Press* the SmartDial™ to enter the setting.
- Turn* the SmartDial™ to adjust the setting.
- Press* the SmartDial™ to accept the change.

3 The Display

The Display is where information can be viewed. The Home Screen is the default display where each point of the Menu System can be viewed. Once a selection has been made within the Menu System, the Display will dim and return to the Home Screen after 30 seconds.

4 Menu System

The 16 points of the clock-face on the Home Screen also represent each of the available settings.

NOTE: These instructions only include functions available for patient use.

Refer to Section 7 for a detailed description of each setting, and PART C of the Diagram Sheet for illustrative examples.

5 Chamber Lid

- Twist the Chamber Lid anti-clockwise to unlock and remove.
- The Chamber Lid must be correctly locked into place for optimal pressure delivery.
- Once the Chamber Lid is locked into place, the handle can be used to carry the device.

6 ThermoSmart™ Breathing Tube

ThermoSmart™ includes a unique heated breathing tube as well as a humidity algorithm which both work to deliver optimal levels of humidity without the side effects of condensation.

CAUTION: Do not use the ThermoSmart™ Breathing Tube if there is any visible damage that exposes copper wires.

7 Elbow

The Elbow is where the breathing tube connects to the device. The Elbow can be rotated from side to side, depending on your preferred sleeping position.

8 Air Filter

The Air Filter can be found at the rear of the device. To access, first pull out the Filter Grill, then pinch the gauze to remove. The Filter Grill must be pushed back into the cavity to keep the Air Filter in place.

9 USB Port

- The USB Port houses the SmartStick™ which must be pushed completely into the USB Port until it clicks into place, flush with the side of the device.
- To remove, push the end of the SmartStick™ for a spring-release retrieval.

NOTE: Only use the SmartStick™ designed and distributed for use with the F&P ICON™.

10 SmartStick™

- The SmartStick™ is used to record sleep data from the device memory when it is inserted into the USB Port of the device and to store music for use with AlarmTunes™. Refer to Section 7 Point 3 for more information.
- The SmartStick™ can be removed and sent to your healthcare provider to view and adjust settings. Changes will be automatically transferred when the SmartStick™ is next inserted into the USB Port of the device.

NOTE: The SmartStick™ must be pulled all of the way out of the USB Port before pushing back in.

11 Power Cord

Four Multi-region plugs are housed in the Carry-Case for personal convenience when traveling. Fit the appropriate plug to the end of the power cord.

12 F&P ICON™ Support Tools CD

The CD contains SmartStick™ Studio Software for transferring music to the SmartStick™, as well as a demonstration video to help with the set-up, and cleaning of the device.

13 Adjustable Straps

Two straps are included with the device. One or both straps can be attached to the Carry-Case according to personal preference when traveling.

14 Carry-Case

The Carry-case has been specially designed to protect the device during transportation. There is a space provided at the handle end of the Carry-Case to store your mask when travelling. Simply remove the Headgear and Elbow to fit. Refer to Diagram 14 in PART B for the correct location of each item when packing.

NOTE: The Carry-Case will not close until the Chamber Lid is removed from the device. This is to ensure water is emptied from the Water Chamber before traveling. Place the Chamber Lid in the partition next to the hinge of the Carry-Case.

15 Use and Care Manual

Read this manual carefully and save these instructions for future reference. If you have any questions please contact your healthcare provider.

Replacement Parts

900ICON108 F&P ICON™ Support Tools CD	900ICON212 Air Filters
900ICON200 Water Chamber	900ICON214 Auto/Premo Lid
900ICON202 SmartStick™	900ICON216 Novo Lid
900ICON204 Elbow	900ICON218 Filter Grill
900ICON206 Outlet Seal	900ICON302 Multi-Region Plug Pack
900ICON208 ThermoSmart™ Breathing Tube*	900ICON306 Carry-Case Straps
900ICON210 Breathing Tube*	900ICON308 Carry-Case

* Applied Parts - to fit .86 inches (22 mm) Conical Connector

7. MENU SYSTEM

Refer to **Part C** of the Diagram Sheet at the beginning of this manual for illustrative examples.

1 On/Off and Ramp

- Press the SmartDial™ to start therapy.
- To activate Ramp, press the SmartDial™ again and hold for 3 seconds. The Ramp symbol will appear .
- To stop and start pressure delivery, press the SmartDial™ while in the Home Screen.

2 Humidity

- Turn the SmartDial™ to the Humidity setting . Press to enter and turn the SmartDial™ again to select the desired humidity level from 0 to 7 (low to high). Press again to accept the change.

NOTE: The default setting is Humidity level 4.

- For additional humidity control, press and hold the SmartDial™ for 3 seconds while in the Humidity setting  and the Boost symbols will appear . Turn to select low, medium or high and press again to accept the change.

NOTE: The default setting is for medium boost humidity.

3 Alarm On/Off

- Turn the SmartDial™ to the Alarm setting , press to enter, then turn to select either On  or Off . Press again to accept the change.
- If Alarm On  is chosen, turn the SmartDial™ to select either the default bell  or the music symbol for AlarmTunes™ , and press again to accept the change.

NOTE: The AlarmTunes™ symbol is only visible when the SmartStick™ is inserted.

- Turn the SmartDial™ to select the desired music track (from 1 to 5 depending on how many tracks are loaded onto the SmartStick™), or select Auto for shuffle and press again to accept the change. Refer to the F&P ICON™ Support Tools CD for SmartStick™ Studio software, as well as instructions for installing and transferring music to your SmartStick™.
- Press the SmartDial™ once to snooze the Alarm for 10 minutes or twice to turn the device and Alarm off.

4 Alarm Time

- Turn the SmartDial™ to the Alarm Time setting , and press to enter. Turn the SmartDial™ to select the hour, press to accept. Turn to select the minutes and press again to accept the change.

5 Alarm Volume

- To adjust the volume of the alarm, turn the SmartDial™ to the Alarm Volume setting  and press to enter. Turn the SmartDial™ to adjust the volume from quiet to loud and press again to accept the change.
- To adjust the click volume of the SmartDial™, press and hold for 3 seconds while in the Alarm Volume setting. The volume symbol  will appear, turn to select the desired volume and press again to accept the change.

6 Clock Time

- Turn the SmartDial™ to the Clock Time setting , and press to enter. Turn the SmartDial™ to select the hour, minutes, 12-hr or 24-hr time. After each selection, press again to accept the change.

7 Brightness

- Turn the SmartDial™ to the Brightness setting , and press to enter. Turn the SmartDial™ to select the Home Screen brightness level and press to accept the change. To adjust the brightness of the display while accessing the Menu System, press and hold the SmartDial™ for 3 seconds and the Active Brightness symbol will appear .
- Turn to select the desired brightness and press again to accept the change.

8 Sleep Data

Turn the SmartDial™ to the Sleep Data setting  and press to enter.

NOTE: The model name will display while in the Sleep Data setting.

Turn the SmartDial™ to scroll through each of the following statistics:

- Total Compliant Hours
- Compliant hours during the last 30 nights
- Percentage Compliant hours over the last 30 nights
- Day that 70% compliance was achieved in 30 consecutive days within the first 90 days of use.

NOTE: A ☺ symbol will appear on the Display Screen when the Sleep Data is within the compliance requirements.

- Set Pressure (for fixed pressure models only)
- Check Sum (⚙ this symbol will only appear if Additional Sleep Data is restricted)
- Exit the Sleep Data setting (⏪ this symbol will only appear if Additional Sleep Data is restricted)

Additional Sleep Data (NOTE: this information may be restricted by your healthcare provider.)

- Leak for last night
- AHI for last night
- 90th percentile pressure for last night
- SensAwake™
- Check Sum
- Exit the Sleep Data setting.

9 SensAwake™ On/Off or Altitude Adjustment

- SensAwake™ is activated as a default setting. To deactivate from models with this feature, turn the SmartDial™ to the SensAwake™ setting  and press to enter. Turn the SmartDial™ to select On ✓ or Off ✗ and press again to accept the change.
- For models without automatic Altitude Adjustment, turn the SmartDial™ to the Altitude setting  and press to enter. Turn the SmartDial™ to select the Altitude level from 0 to 10,000 ft and press to accept the change.
- To choose meters, press the SmartDial™ and hold for 3 seconds while in the Altitude Adjustment setting. When the symbol starts to flash, turn to select meters and press to accept the change. Press the SmartDial™ again to enter the menu, and turn to adjust to the desired Altitude level. Press again to accept the change.

10 ECO Mode

ECO Mode allows for the device to operate on a reduced power (75 W).

NOTE: Humidity output will be significantly reduced in this mode.

Situations for use in this mode might include long-haul air flights and camping.

- Turn the SmartDial™ to the Eco setting **ECO** and press to enter. Turn to select On ✓ or Off ✗ and press again to accept the change.
- The Eco symbol **ECO** will be displayed continuously on the Home Screen while the Economy mode is activated.

11 Display Pressure

The Display Pressure setting provides two display options for the Home Screen. Time is the default setting, however this can be changed to show pressure (cmH₂O or hPa) if required.

- Turn the SmartDial™ to the Display Pressure setting  and press to enter. Turn to select On ✓ or Off ✗, and press again to accept the change.
- If a different unit of measure is required, press the SmartDial™ and hold for 3 seconds until the symbol  appears. Turn the SmartDial™ to select either cmH₂O or hPa and press again to accept the change.

12 Service Setting

The Service setting will only appear if an error has occurred. If the Error symbol **Error** flashes on the Home Screen, turn the SmartDial™ to the Service setting  and write down the number displayed. Contact your healthcare provider for further instructions.

8. COMMUNICATING SLEEP DATA

Telephone Reporting

Many healthcare providers will ask for sleep data over the telephone to track CPAP therapy progress. Follow the steps outlined under Sleep Data (Section 7, Point 8) which explains how to access this data.

Downloading Data to the SmartStick™

If your healthcare provider requests a copy of your sleep data from the SmartStick™ please follow the steps below to download this information:

- (1) If the SmartStick™ is already inserted into the USB Port of the device, check that the pressure is off and there is no SmartStick™ symbol  flashing on the Home Screen. The SmartStick™ can now be safely removed by pushing the end for a spring-release retrieval.
- (2) If the SmartStick™ has been removed, ensure the pressure is off before inserting the SmartStick™ back into the USB Port of the device. Wait until the SmartStick™ symbol  has stopped flashing on the Home Screen before removing the SmartStick™ again.
- (3) Data from the SmartStick™ should then be sent in the manner requested by your healthcare provider.

NOTE: The flashing symbol indicates that data is being transferred from the device memory to the SmartStick™. Removal of the SmartStick™ during this time may cause data corruption.

9. CLEANING AND MAINTENANCE

The device should be cleaned as follows:

- Unplug the device from the power supply.
- Wipe the exterior and basin of the device (when the Water Chamber has been removed) with a clean, damp (not wet) cloth and mild dishwashing detergent. Do not use harsh abrasives or solvents, as these may damage the device.

Daily

Clean the Water Chamber and the breathing tube:

- Disconnect the breathing tube from the mask and the device Elbow.
- Clean the breathing tube with warm soapy water, rinse thoroughly and hang with both ends pointing to the floor to dry.
- Remove the Water Chamber and discard any remaining water.
- The Water Chamber can be washed in a domestic dishwashing machine, or cleaned and washed with mild dishwashing detergent, then rinse and dry.

Weekly

Thoroughly clean the Water Chamber:

- Soak the inside of the Water Chamber for 10 minutes in a solution of one part white vinegar to two parts water.
- Empty the vinegar solution and rinse the Water Chamber thoroughly with water.

Replace the Air Filter when it becomes significantly discolored, at least once every three months or after 1,000 hours' worth of running time.

With a regular cleaning routine, it is recommended to replace your breathing tube and Water Chamber every six months. However equipment should be replaced immediately if there is any sign of cracking, discoloration or leaking.

NOTE: This device is not serviceable and should not be disassembled. Refer all therapy and equipment queries to your healthcare provider.

10. PRODUCT SPECIFICATIONS

DIMENSIONS:	6.3 H x 6.7 W x 8.7 D in. (160 H x 170 W x 220 D mm)
WEIGHT:	4.8 lb (2.2 kg) Packed weight 8.7 lb (4.0 kg)
PERFORMANCE:	
Pressure Range:	4 to 20 cmH ₂ O (in the unlikely event of fault conditions, pressure may reach up to 26 cmH ₂ O)

Maximum Flow Rates					
CPAP Pressure Setting (cmH ₂ O)	4	8	12	16	20
Measured pressure at patient connection port (cmH ₂ O) - AUTO/PREMO	62	102	129	149	143
Measured pressure at patient connection port (cmH ₂ O) - NOVO	47	52	57	61	65

Dynamic Pressure Stability					
CPAP Pressure Setting (cmH ₂ O)	4	8	12	16	20
Pressure difference (cmH ₂ O) at mask connection port - AUTO/PREMO	0.6	0.6	0.6	0.7	0.8
Pressure difference (cmH ₂ O) at mask connection port - NOVO	1.2	1.2	1.3	1.5	1.5

Static Pressure Stability		
	AUTO/PREMO	NOVO
Change at connection port at a pressure setting of 10 cmH ₂ O	0.2	0.4

HUMIDITY: **Maximum humidity** = 35 mg/L (BTPS), 85% RH at 10 cmH₂O, with humidity level 7 and Boost level high.

Typical humidity = 21 mg/L (BTPS), 78% RH at 10 cmH₂O, with humidity level 4 and Boost level medium.

GAS TEMPERATURES: Maximum = 100 °F (38 °C)

NOISE LEVEL: Sound Pressure Level < 29 dBA
Average Sound Power Level < 37 dBA
420 ml up to the maximum fill line

WATER CHAMBER VOLUME:

STANDARDS COMPLIANCE: IEC60601-1: 1988 + A1 - A2
AS/NZS 3200.1.0: 1998
EN60601-1: 1990

DATA RECORDING: The SmartStick™ records up to 365 days of summary efficacy data, 30 days of detailed efficacy data (AHI, Leak), 7 days of high frequency Flow, Pressure and Leak data (50 Hz), as well as cumulative data from first use.

Electrical Ratings		
Rated supply voltage	Rated current input	Rated supply frequency
110 - 115 V	1.27 A (1.43 A Max)	50 - 60 Hz/400 Hz
220 - 240 V	1.07 A (1.21 A Max)	50 - 60 Hz

STORAGE AND TRANSPORT: The device should always be stored and transported within the following temperatures:

Transport: 14 to 140 °F (-10 °C to 60 °C)

Storage: 14 to 140 °F (-10 °C to 60 °C)

CAUTION: Ensure water is emptied from the Water Chamber before traveling. The Carry-Case will not close until the Chamber Lid is removed from the device.

OPERATING CONDITIONS: **Ambient Temperature:** 41 to 95 °F (5 to 35 °C)

Humidity: 5 to 95% Relative Humidity

Altitude: 0 to 9,000 ft (0 to 3,000 m)

NOTE: Above 4500 ft (1500 m) the maximum operating pressure will be reduced.

11. DEVICE AND CONSUMABLES DISPOSAL INSTRUCTIONS

Unit Disposal Instructions

This unit contains electronics. Please do not discard as regular rubbish. Dispose of electronics according to local guidelines.

Consumable Disposal Instructions

Place the mask, breathing tube and Water Chamber in a rubbish bag at the end of use and discard with normal rubbish.

12. TROUBLESHOOTING

If you feel that your device is not operating correctly, please refer to the following suggestions. If the problem persists, please consult your healthcare provider. Do not attempt to repair the device yourself.

Problem	Possible Cause	Solution
There is no visible display or power	The plug may not be correctly attached to the power cord	Ensure the correct plug from the Carry-Case is attached to the power cord, and is connected to the power supply.
Insufficient air is delivered from the device	The Chamber Lid may not be correctly fitted to the device	Refer to Section 5 and Part B for detailed instructions for fitting the Chamber Lid.
"Error" is displayed on the Home Screen of the device	An error may have been detected with the device	Turn the SmartDial™ to the Service setting  and press to enter. Write down the number displayed and contact your healthcare provider for further instruction.
The Carry-Case won't close	The Chamber Lid may not be removed from the device	The Carry-Case has been designed to protect the device during travel. The Chamber Lid must be removed before closing the Carry-Case, and water must be removed from the Water Chamber.
AlarmTunes™ does not work	The Alarm may be set to "off" The Alarm may be set to "buzzer" The SmartStick™ may not have been inserted into the USB Port correctly	Refer to Section 7, point 3 for instructions on the Alarm setting Once SmartStick™ Studio software has been installed onto your computer and your favorite songs selected, download onto your SmartStick™ and then ensure the SmartStick™ is inserted into the USB Port of the device.

13. F&P ICON™ MODELS AND FEATURES MATRIX

The model name of each device is located on the base, as well as in the Sleep Data setting of the Menu System. Refer to Section 7, Point 8 for more information.

Features	AUTO	PREMO	NOVO
Fully integrated CPAP with humidification	•	•	•
ThermoSmart™ Heated Breathing Tube	•	•	•
Auto Adjusting Pressure	•		
*Efficacy Reporting (AHI, Leak)	•	•	
Compliance Reporting	•	•	•
SmartStick™ Removable Media	•	•	•
SensAwake™ Pressure Relief	•		
Proportional Ramp	•	•	•
Auto-Altitude Adjusting	•	•	Manual
Leak Compensation	•	•	

* **NOTE:** The display of efficacy data may be restricted by your healthcare provider.

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