

Voyage™

Go Together. Grow Together!

The universal Zippie[®] Voyage[™] early intervention device combines Baby Jogger[™] stroller expertise, Zippie's incredibly versatile seating for growing babies, and one of the largest selections of functional accessories in the industry!

Baby Jogger City Select Base Zippie teamed up with Baby Jogger to give you an intuitive, easy-to-use stroller option that makes no compromises when it comes to what's best for baby.

Telescoping Stroller Handle Responsive, one-hand steering that adjusts to caregiver heights.

Canopy with Mesh Window When the seat is facing forward, open the mesh window to keep an eye on your little one.





Cross Over

Center Stripe



Quick-Access, Hand-Operated Wheel Locks When juggling kids, bags, or medical equipment, locking the wheels is as simple as flipping one lever on the stroller handle.

Functional Accessories

The Voyage is available with options and accessories designed to make daily tasks a breeze, so you can delight in the little joys of caring for your family.



WC-19 Transit Approved Approved for use as an occupant's seat in a motor vehicle when used with an approved Wheelchair Tiedown and Occupant Restraint System in accordance with SAE (J2249).

Folding Booster Base

(Height Adjustable

O2 Cylinder Holder

Carry-All Storage Basket Conveniently store up to 15 lbs. of your essentials.

Patented Quick-Fold Technology Folds into a compact package in one simple step so you can be on your way.

Puncture-Resistant Rear Wheels Travel confidently knowing the Voyage is up for the trip.

Swivel Front Wheels with Suspension and Lock-Out Glide over and around obstacles while keeping baby calm, secure, and comfortable. Lock-out the front wheels for smooth tracking during long-distance strolls.







Battery and Vent /

Utility Under

Frame Trav



Utility Hook Upper Extremity Support Tray (Angle and Depth adjus

Dovetail Vertical Frame

Mounted Hardware for

LTV/Trilogy Ventilators

Rear Add On Accessory Platform Seating for Babies on the Grow! When it comes to seating, the Voyage team's number one priority is positioning your child for healthy growth and development from babyhood through their toddler years.



With Hip Blocks Without Hip Blocks

Standard Planar Seat

- Contoured foam
- · Fixed Ischial Ledge
- Built-in Depth Growth



Adjustable Contour Seat

- Adjustable Pre-Ischial Ledge
- Lateral Hip Supports
- 3" Seat Depth Growth Kit
- Removable neoprene cover

Height and Angle Adjustable Footrest with optional foot straps and shoe holders



Lateral Trunk Support

Positioning Accessories



Butterfly Harness

Positioning Belt

Tilt/Recline

to assist with feeding, digestion, respiratory function, and visual orientation.

Forward (Anterior) Tilt Assists with positioning to encourage achievement of developmental milestones such as bringing hands to the mouth, reaching for objects, or engaging in activities.



- T Back
- Short (7''– 13'')
- Tall (10"-17")
- 1" or 2.5" thick foam

Whitmyer Head Supports





Plush

Medial Support

Contour Cradle

Adjust-a-Plush

Foot Straps



Improves positioning



Easy Access Tilt/Recline Release Tilting and reclining the seat is as simple as activating an easy-to-reach trigger handle behind the backrest.

Reversible Seat

Face the seat toward vou to interact with your child, or forward so your child can explore their surroundings.



Color Options



Specifications

| TECHNICAL SPECIFICATIONS | |
|--|----------------------------------|
| Frame Style: | Early Intervention Device |
| HCPCS Code: | E1232 E1234 (without seating) |
| Chair Weight: | Base only: 23 lbs. |
| Frame Material: | Aluminum |
| User Weight Capacity: | 60 lbs. |
| Transit Approved: | Yes (WC-19) |
| | |
| FRAME AND FEATURE MEASUREMENTS (Dependent on Configuration) | ; |
| | 24" |
| (Dependent on Configuration) | |
| (Dependent on Configuration) Overall Width: | 24" |
| (Dependent on Configuration) Overall Width: Seat Widths: | 24" 6" - 13" |
| (Dependent on Configuration) Overall Width: Seat Widths: Seat Depths: | 24" 6" - 13" 5" - 13" |

| OPTIONS & ACCESSORIES | |
|-----------------------|--|
| Frame Type: | Til l Tilt & Recline |
| Caster Wheels: | 7" |
| Rear Wheels: | 12" Pneumatic with Airless Inserts |
| Backrest Type: | T-Back Cushion 1" Thick T-Back Cushion 2.5" Thick |
| Seating: | Standard Planar Seat without Hip Blocks Standard Planar Seat with Hip Blocks Contoured Wedge Seat Contoured Anti-Thrust Seat Adjustable Contour Seat |

The Zippie Voyage

How to Unfold & Fold the Zippie Voyage



How to Tilt, Recline, & Reverse the Zippie Voyage Seat



QR Codes to Voyage Videos



Customer Service: 800.333.4000 www.SunriseMedical.com