

SafeFree® - The Patented Side Guard Concept from Wissner-Bosserhoff

No protection needed

Barrier-free; no middle post and no central gap protector



In lowered position, the divided side guards enable unhindered access to the resident for daily nursing tasks, changing wound dressings or bed sheets and linen. No need to take off and store interfering accessories like a middle post or fixed side guards.

Soft protection – no FDM (freedom-depriving measures)

Orientation instead of freedom-depriving measures



For safety at night, it is often sufficient to combine the low position of the bed with the first height level of the head-end side guard. This “soft protection” avoids unnecessary barriers and reduces the risk of falls sustainably. The side guards provide the resident with orientation without restriction.

3/4 protection – no FDM (freedom-depriving measures)

Mobilization with adaptation to resident height



The 3/4 protection creates a deliberate bed egress gap at the foot end for the reduction of freedom-depriving measures (FDM), while at the same time maintaining a high safety standard. By adjusting the bed and side guard height, sit height can be adapted to residents of every size and height.

Full protection – for mattress heights from 4.7” to 7”

Safety



The third height level is used to comprehensively secure the resident in bed. In the top position, the two parts of the side guard close together so tightly that they are just as safe as a continuous side guard. Stocking up on accessories for closing a central gap when using an abdominal belt for restraint is not necessary.