



PRESSURE
ULCER
TREATMENT

EQUAL PRESSURE DISTRIBUTION AND PERFECT SUPPORT FOR THE ASYMMETRIC BODY

Pre-ischial bar prevents
sliding forward | Converts
shear into normal forces

Thigh support

Total immersion for equal pressure
distribution | Support IT-coccyx area |
Managing pelvic asymmetry.



An individual solution for asymmetry, even the extreme kinds

Total immersion for equal pressure distribution

Prevents sliding while minimizing shear

The Adjuster is extremely suitable for use with asymmetric body contours and amputations. It automatically manages small obliquities and is individually adjustable to support more extreme asymmetry. When you are looking for a high level of comfort, optimal skin protection and stability: the Vicair Adjuster is always a safe choice. The design of this wheelchair cushion allows total immersion, to redistribute pressure equally. At the same time it prevents you from sliding forward. As a result the Adjuster offers you optimal protection against skin breakdown.

PRODUCT PERFORMANCE

Skin protection / Anti decubitus	▲▲▲▲
Pressure redistribution	▲▲▲▲
Shear reduction	▲▲▲▲
Micro climate control	▲▲▲△△
Impact & vibration damping	▲▲▲▲

Positioning

Lateral stability	▲▲▲▲△
Frontal stability	▲▲▲▲
Sliding prevention	▲▲▲▲
Allows manual adjustments	▲▲▲▲
Allows easy transfers	▲▲▲▲△

Maintenance & Safety

Durability	▲▲▲▲
Reliability	▲▲▲▲
Functionality when misused	▲▲▲△
Low maintenance level	▲▲▲△
Level of hygiene	▲▲▲△

Functionality & Application

Comfort	▲▲▲▲
Low cushion weight	▲▲▲▲
Available size range	▲▲▲▲
Maximum user weight	▲▲▲▲
Low seat to floor height	▲▲▲▲
Auto adjustments	▲▲▲△

USED WITH

Hemiplegia (CVA), amputations, foot propulsion, kyphosis, pelvic obliquity/tilt/rotation, Spina Bifida, spasms, Fibrodysplasia Ossificans Progressiva (FOP)

TECHNICAL SPECIFICATIONS

Height	6 and 10 cm
Average weight	534-769 grams
Max. user weight	250 kg

CUSHION SIZES

Standard	
On request	
Inch	14 15 16 17½ 18 19 20 21 22 23 24
CM	35 38 40 43 45 48 50 53 55 58 60
14	35
15	38
16	40
17½	43
18	45
19	48
20	50
21	53
22	55
23	58
24	60

en-rev 04/2016