


# Dimensions of the Wave

User dimensions (inches)	Z210 small	Z220 medium	Z230 large
Height	up to 46	40–56	50–74
 <b>Key user dimension: height</b> Select the appropriate bath chair by the user's overall height. Choose the size that allows for growth. <b>Important:</b> Make sure that seat depth, width and back height are adequate. The user's weight must not exceed the maximum working load.			
Item dimensions (inches)	Z210 small	Z220 medium	Z230 large
Seat depth	10½	13	17
Seat width	15	15	17½
Back height	23	26½	32
Seat height	2–7	2–7	2–7
Seat height (incl. tub stand)	12½–18½	12½–18½	12½–18½
Seat height (incl. shower stand)	27½	27½	27½
Seat height (incl. tub transfer base) with long leg inserts	21–26 24½–29½	21–26 24½–29½	21–26 24½–29½
Length of mesh area when flat (without calf rest)	35	41	50½
Length of mesh area when flat (with calf rest)	42½	49	61
Width of base	18	18	21
Length of base when flat	29½	29½	29½
Tub stand height (legs flat)	8	8	8
Tub stand length (legs flat)	39	39	39
Tub stand width (legs flat)	17½	17½	17½
Shower stand height	22	22	22
Shower stand width	27	27	27
Shower stand length	34	34	34
Tub transfer base footprint	28 x 30	28 x 30	28 x 30
Max tub height for transfer base with long leg inserts	20 24	20 24	20 24
Maximum working load (lbs)	75	150	250
Back and seat angle adjustments		0°, 15°, 30°, 45°, 60°, 75°	



**Without calf rest**  
(Z269)



**With calf rest**  
(Z267)

## Conversion kits

are the smart solution for growing children. For a fraction of the cost you can convert a small to a medium bath chair.

We'll send you a new larger fabric cover, back and seat frames and chest strap with instructions on how to remove the old and install the new parts on your existing base.



The Wave has a washable fabric cover that is easily removed for cleaning.



All three sizes are available in two beautiful colors, blue and pink.