

ULTIMATE FREEDOM FOR THE MOBILE USER



Stays with you during transfers

Skin protection

Extremely light to support optimal mobility

Active as you are, not all of your activities require or allow the use of a wheelchair. The Vicair AllRounder keeps you comfortable and protects your skin. As it is easily buckled onto your body, it follows you around wherever you go.

Use it during your (wheelchair) sports-activities, to sit in a normal chair or to play on the ground. Experience the freedom!

PRODUCT PERFORMANCE

| | |
|----------------------------------|-------|
| Skin protection / Anti decubitus | ▲▲▲△△ |
| Pressure redistribution | △△△△△ |
| Shear reduction | △△△△△ |
| Micro climate control | △△△△△ |
| Impact & vibration damping | ▲▲▲△△ |

| | |
|---------------------------|-------|
| Positioning | △△△△△ |
| Lateral stability | △△△△△ |
| Frontal stability | △△△△△ |
| Sliding prevention | △△△△△ |
| Allows manual adjustments | △△△△△ |
| Allows easy transfers | △△△△△ |

| | |
|----------------------------|-------|
| Maintenance & Safety | ▲▲▲△△ |
| Durability | ▲▲▲▲△ |
| Reliability | ▲▲▲▲△ |
| Functionality when misused | △△△△△ |
| Low maintenance level | ▲▲▲▲△ |
| Level of hygiene | ▲▲▲▲△ |

| | |
|-----------------------------|-------|
| Functionality & Application | ▲▲▲△△ |
| Comfort | ▲▲▲▲△ |
| Low cushion weight | ▲▲▲▲△ |
| Available size range | ▲▲▲△△ |
| Maximum user weight | ▲▲▲▲△ |
| Low seat to floor height | △△△△△ |
| Auto adjustments | △△△△△ |

USED WITH

Regular transfers on and off the wheelchair | When seated on hard surfaces for long periods in time.

TECHNICAL SPECIFICATIONS

| | |
|------------------|-----------|
| Height | n.a. |
| Average weight | 800 grams |
| Max. user weight | 100 kg |

CUSHION SIZE

Available in one size