

Invacare® Pronto® M41™ Power Wheelchair



Yes, you can.®

Pronto M41™





Yes, you can.®

The Invacare® Pronto® M41™ power wheelchair is the latest offering in the sleek and versatile Pronto® series of power wheelchairs. The M41™ features the latest technology in electronics and in-line motors. This combination provides increased efficiency and battery range, as well as an increased top speed. The M41 also features a center-wheel drive platform, simplified SureStep®, and six wheels on the ground providing stability and maneuverability similar to other Pronto® models. All of these features are packed into a simple, clean-looking, and easy-to-service compact frame that provides mobility with style for people on the go!



Pronto® M41™ with Fold-Down Seat



Pronto® M41™ with Semi-Recline Office-Style Seat



Pronto® M41™ with Solid Seat Pan.

Features


- Compact size for a tight turning radius
- High efficiency in-line motors for greater speed and longer battery range
- SureStep® suspension with six wheels on the ground for smooth driving and optimum stability
- Charger connection conveniently located on joystick
- Stylish design
- 3-amp off-board charger
- Three seat options
 - ✓ Medium-back fold-down seat
 - ✓ High-back semi-recline office-style seat
 - ✓ High-back semi-recline van seat with solid seat pan
- Auto style positioning strap
- U1 batteries
- 300 lb. weight capacity

Colors Available

Candy Red



Deep Blue



Pronto® M41™ Specifications

Overall dimensions (measured with fold-down seat)			
Length	34"		
Width (including joystick)	24"		
Height	33"		
Seat-to-floor	19"		
Seat dimensions	(fold-down seat)	(semi-recline seat)	(semi-recline seat w/solid seat pan)
Width	18"	16" - 20"	16" - 20"
Depth	17"	16" - 18"	16" - 20"
Seat angle	92°	80° - 115°	90° - 115°
Weight capacity	300 lb.		
Wheel sizes			
Front/rear casters	6" x 2"		
Drive wheels (Flat-free)	10.5" x 3.5"		
Maximum speed	5 mph		
Motor	250 watts		
Incline capability	7° grade		
Turning radius			
Front	17.75"		
Rear	19.5"		
Weight	(fold-down seat)	(semi-recline seat)	(semi-recline seat w/solid seat pan)
with U1 batteries	161 lb.	175 lb.	180 lb.
w/out batteries and foot plate	110 lb.	124 lb.	129 lb.
Ground clearance	2"		
Battery	U1		
Battery charger	3 amp off-board		
Per-charge battery range	Up to 17 miles*		
Limited warranties			
Base frame	60 months		
Seat frame	36 months		
Electronics	12 months		
Motor/gearbox	12 months		

* Travel range will vary depending on user weight, speed, battery charge and terrain.

Do not operate on roads, streets or highways.

Invacare Corporation
www.invacare.com

USA

One Invacare Way
Elyria, Ohio
44035-4190
(800) 333-6900

Canada

570 Matheson Blvd. E.,
Unit 8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324

©2008

Invacare Corporation.
All rights reserved.
Trademarks are identified by the symbols ™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted.
Specifications are subject to change without notification.

Form No. 06-033
Rev. 10/08