



Powerstep®



Covered under
US Patent
D661,076



PINNACLE MAXX FULL LENGTH

MATERIALS

- Polypropylene support shell with 2 degree medial post and deeper heel cup
- EVA base, Charcoal Grey
- Polyester fabric top cover with anti-microbial protection, Maroon
- VCT® top layer, Charcoal Grey, 1/8"

FEATURES & BENEFITS

- Strong, medical-grade foot support calibrated with flexibility
- Deeper heel cup and 2-degree medial post for improved foot control and positioning
- Double layer, shock absorbing EVA/VCT® cushion casing for enhanced comfort
- Friction reducing anti-microbial top fabric

INDICATIONS: Over-Pronation (Moderate to Severe); Plantar Fasciitis/Heel Pain; Foot, Arch and Heel Pain; Calluses; Hammertoes and Claw Toe; Bunions; Neuromas; Ankle and Knee Pain

IDEAL SHOE TYPE: Roomier Dress and Casual Shoes, Athletic Shoes

Product Code	Size	Men's Size	Women's Size	UPC Code
5015-01Asm	Asm	3 - 3 ½	5 - 5 ½	749852100161
5015-01A	A	4 - 4 ½	6 - 6 ½	749852100260
5015-01B	B	5 - 5 ½	7 - 7 ½	749852100369
5015-01C	C	6 - 6 ½	8 - 8 ½	749852100468
5015-01D	D	7 - 7 ½	9 - 9 ½	749852100567
5015-01E	E	8 - 8 ½	10 - 10 ½	749852100666
5015-01F	F	9 - 9 ½	11 - 11 ½	749852100765
5015-01G	G	10 - 10 ½	12	749852100864
5015-01H	H	11 - 11 ½		749852100963
5015-01I	I/J	12 - 13		749852110160
5015-01K	K	14 - 15		749852110269
5015-01L	L	16+		749852110252

WEARING/FITTING DIRECTIONS: Remove the shoe's existing insoles if possible and replace with your Powerstep® insoles. Trim edges if needed to fit. If necessary, wear a few hours a day until feet are accustomed to better support, and then wear full time. Suggest wearing Powerstep® insoles in all your shoes, even slippers, to give continuous foot support.

CARE INSTRUCTIONS: Spot clean with a mild detergent. Do not immerse in water. Air dry only. Replace your Powerstep® insoles as needed depending on wear.

WARNING: Not intended for persons with severe foot deformities or circulatory disorders resulting in infections or foot ulcerations. Consult your podiatrist or medical professional before using. Always consult a medical professional when experiencing persistent pain, discomfort or swelling.



888.237.3668 · 513.825.1888
Email: info@powersteps.com
www.powerstep.com
Stable Step, LLC, West Chester, OH
Made in the USA with US and imported materials.

Powerstep® is a registered trademark of Stable Step, LLC.
VCT® is a registered trademark of Remington Products, Inc.

HCPCS L3040 Approved