

EasyStand Bantam medium



SIT + STAND + SUPINE
STAND BIGGER KIDS

The Bantam

Hip-Knee Flexion Transfers

- Allows standing pivot transfers, lateral seat to seat, supine flat-to-load, supine with hip and knee flexion transfers, or any position in between
- Compatible with patient lifts
- Each child is different and the Bantam provides choices when transferring

Gradual ROM with Hip-Knee Flexion

- Can start in a hip-knee flexion position and gradually increase lower extremity extension
- By maintaining a sustained, weight-bearing stretch throughout the standing session, range of motion (ROM) can increase over time for many children



Full Support with Sit to Stand

- Full body support system with foot plates, knee pads (anterior), seat (posterior) and chest pad (anterior) that supports the individual in all positions from sitting, supine, and standing
- Children appear to adapt more quickly to standing and weight-bearing when using a sit to stand device due to the extra support and sense of security

Accommodating Contractures

- Adjustability creates the opportunity to stand children with positioning challenges
- No special modifications are needed to accommodate children with contractures and leg length discrepancies
- Match the hip and knee flexion of the child and then adjust as range of motion in the hips and knees improves



Advantage

Upright Standing

- Upright standing is a natural position with the head, hips, knees and ankles all in vertical alignment
- When a child stands upright, they are likely to feel more balanced and are in an optimal position to allow strengthening of the neck muscles
- The option to easily go from upright to supine encourages longer periods of standing

Custom Fit for Better Standing

- Modular design works for children who require minimal positioning and has options for those who require maximum positioning and support
- Over 50 positioning and support options that help keep the most challenging cases safe, aligned and secure

Positioning Chair Capabilities


- Can do anterior tilt positioning, posterior tilt in space, and a true full recline similar to positioning/activity chairs
- Offers terrific pelvic support and stability by keeping the child positioned in the seat with knee pads, hip supports, and lap belts
- More position changes with fewer transfers means less lifting and decreased risk of injury for caregivers

Plenty of Growth

- Substantial growth built into the frames with tool free adjustments
- The Bantam medium has up to 18" of height growth
- Seat depth: 15"-20"
- Seat to Foot Plate Range: 7"-18.5"



Options




PY5554
Supine

Allows the Bantam to be used as a supine or sit to stand stander. User can be transferred in flat-to-load, seated, or hip-knee flexion position. Adjusts from 0-90° supine. Head support required. Not available off unit. Shadow tray recommended.




PY5632
Adjustable/Removable Actuator Handle

Removable actuator handle adjusts the distance of reach to handle to meet user's needs. Not available with Pow'r Up Lift.




PY5608
Black Molded Tray for Swing-Away

Provides user with a tray surface and anterior support in the standing position. Includes chest pad, adjustable in depth and height. Tray size is 21"x24" (53x61cm).



PY5618
Oversized Angle Adj. Tray for Swing-Away

Allows the tray to swing away providing more room for front transfers. Tray size is 29"x22" (74x56cm).




PNG30000
Large Contoured Chest Pad

Attaches to tray or no-tray bracket to provide anterior support with a minimal lateral contour. Oval pad is 9.5"Hx15"W (24x38cm). Not available with elbow pads.



PY5600
Roho® Knee Pads

Includes Roho® Dry Floatation® technology. Width range from center to center is 5"-9" (13-23cm). Depth range from front of seat is 3"-8" (8-20cm). Height range is +/- 1.5" (+/- 4cm) from seat pivot.




PY5574, PY5576
Contoured Back 16"H, 20"H

Contoured back provides midline support. Top of back is 16" (41cm) or 20" (51cm) from seat.



PNG30171, PNG30021
Y-Style Chest Vest-Medium, Large

Neoprene®. Medium size is 9.5"Lx9"W (24x23cm). Large size is 11"Lx9.5"W (32x26cm). Back required.



PY5643
Angle Locator

Helps therapists monitor standing progress. Attaches on the back with magnet.



PY5645
Easy-Adjust Seat Depth

Roller bearing seat depth that allows for smoother seat depth adjustment. Not available off unit.



PNG50044
Quad Grip Handle Extension

Attaches to the actuator handle for users with limited dexterity. Adjusts for various hand sizes. Not available with Pow'r Up Lift.



PY5610
Clear Tray for Swing-Away

Provides user with an acrylic tray surface and anterior support in the standing position. Includes chest pad, adjustable in depth and height. Tray size is 21"x24" (53x61cm).



PNG50042
Hand Grips

Provides additional arm and hand support and alignment. Easily movable for precise arm positioning. Attaches with suction cups. A clear acrylic tray is required.



PNG30295
Padded Tray Cover

A padded tray cover that attaches to the black molded, clear, angle adjustable, or shadow tray. Not available with oversized clear tray or elbow pads.



PY5596
Swing-Away Knee Pads

Push of a button swings knee pad away providing more room for transfers. Multi-adjustable knee pads are adjustable in height +/- 3" (+/- 8cm), width 7"-9" (13-23cm) in addition to depth range 3"-6" (8-15cm).



PY3008
Velcro® Positioning Belt

2" (5cm) wide positioning belt has a D-ring Velcro® closure, provides hip stability for user. Fits hip circumference up to 66" (168cm). Back required.



PY5562
Lateral Supports

Provides midline positioning for users with limited trunk control. Width adjustment range is 9"-14.5" (23-37cm). Pad size is 4"Hx6"W (10x15cm). Not available with contoured back.



The Bantam medium base model comes with two 5" rear locking casters, two 5" front wheels, manual hydraulic actuator with handle, planar seat, knee pads, height-only adjustable foot plates, no-tray chest pad, black upholstery, swing-away front (supine not included).



PY5568
Front Swivel Casters for Shadow Tray

Allows for easier movement of stander. Replaces front wheels. Extends the footprint of the unit to 28.25"x47" (72x119cm).



PNG50070
T-Style Handle Extension

Attaches to the actuator handle for users with limited range of motion. Makes the first few pumps easier and descending more accessible. Adds additional 7" (18cm) in each direction. Not available with Pow'r Up Lift.



PY5612
Oversized Clear Tray for Swing-Away

Large acrylic tray allows more room for reading and functional activities in standing position. Tray size is 29"x22" (74x56cm).



PT50156
Elbow Pad-9.25" Cutout

Additional padded support surface on the black molded or clear tray. Not available for oversized clear tray. 9.25" cutout width. Overall width 21.5" (55cm), pad extends back 3.75" (10cm) posterior from current table edge.



PY5586
Multi-Adjustable Foot Plates

Foot plates are adjustable forward/aft +/- .75" (2cm), plantar/dorsi +/- 20°, toe-in/-toe-out, swivel 30° left/right, and height.



PY5570, PY5572
Contoured Seats-Narrow, Wide

Contoured positioning seat features form-to-fit upholstery with built in hip support. Narrow seat is 11"-13"W (28-33cm). Wide seat is 12"-14"W (30-36cm).



PY3014
Positioning Belt w/Airline Style Buckle

2" (5cm) wide positioning belt with airline style buckle provides hip stability. Fits hip circumference up to 66" (168cm).



PY5626, PY5628, PY5630
Head Support-Extra Small, Small, Large

Height range from seat with 16" planar back is 18"-27" (46-69cm), 20" planar back is 22"-31" (56-79cm). Depth range is -2.5" to +1.5" (-6 to -4cm) from backrest.



PT50170
Front Swivel Casters for Swing-Away Front

Allows for easier movement of stander. Replaces front wheels. Extends the footprint of the unit to 29.5"x47.5" (72x121cm).



PY5604
Black Molded Shadow Tray

Shadow tray is accessible in all positions and provides anterior support as user moves from sitting to standing. Adjustable in depth, height, and tray angle. Tray size is 21"x24" (53x61cm).



PY5614
Black Molded Angle Adj. Tray for Swing-Away

Adjusts from 0-35° to provide forearm support and angle activities closer to the user in the standing position. Tray size is 21"x24" (53x61cm).




PT50160
Elbow Pad-11" Cutout

Additional padded support surface on the black molded or clear tray. Not available for oversized clear tray. 11" (28cm) cutout width. Overall width 21.5" (55cm), pad extends back 3.75" (10cm) posterior from current table edge.



PY5636
Foot Straps

Velcro® and D-ring adjustment holds feet in alignment with multiple attachment slots for proper foot positioning.



PY5560
Hip Supports

Supports the hips in symmetrical alignment from sitting to standing. Adjustable width from 9"-15" (23-38cm) with a 1" (2.5cm) height range. Pad size is 4"x6" (10x15cm). Not available with contoured seat.



PY5620
Chest Strap

Velcro® and D-ring adjustment provides moderate support, back required. Circumference up to 42" (107cm). Pad size 5"Wx8"L (13x20cm).



PNG50040
Push Handles

Makes pushing the stander easier for the caregiver. Back required.



PY5564
Pow'r Up Lift

Brings the user to the standing position with the touch of a button. Includes light touch controller, rechargeable battery and charger, Linak® actuator and an emergency release.



PY5606
Clear Shadow Tray

Shadow tray is accessible in all positions and provides anterior support as child moves from sitting to standing. Adjustable in depth, height, and tray angle. Tray size is 21"x24" (53x61cm).



PY5616
Clear Angle Adjustable Tray for Swing-Away

Adjusts from 0-35° to provide forearm support and angle activities closer to the user in the standing position. Tray size is 21"x24" (53x61cm).




PY5582
Elbow Pad-13" Cutout

Additional padded support surface on the black molded or clear tray. Not available for oversized clear tray. 13" (33cm) cutout width. Overall width 21.5" (55cm), pad extends back 3.75" (10cm) posterior from current table edge.



PY5638
Secure Foot Straps

Velcro® and D-ring adjustment holds feet securely in place at ankle and toes with multiple attachment slots for proper foot positioning.



PY5530, PY5532
Planar Back 16"H, 20"H

Flat back provides additional support. Top of back is 16" (41cm) or 20" (51cm) from seat.



PNG30170, PNG30056
X-Style Chest Vest-Medium, Large

Neoprene®. Medium size is 9.5"Lx9"W (24x23cm). Large size is 11"Lx9.5"W (32x26cm). Back required.



P80809
Tool Pouch

Holds the assembly tools provided with the stander. Attaches to unit with Velcro® straps. Tools not included.

Promoting Activity in a Bantam

- Keeps a child moving all day long
- Infinite positions can make activities easier and more comfortable
- Being able to move the child into several different positions with one product can create opportunities for better health, cognitive growth, and balance

Feeding Chair Possibilities

- Sitting upright is often not the best position for head control for feeding and partial or supine standing can be the answer
- Provides many position options including sitting, tilt, recline, and standing to find the optimal feeding position



Inclusion in School and Home

- The child can stand when others are standing, sit when they are sitting, and use the stander under the same tables, up next to the white board or with other education tools
- Serves as a positioning, feeding, or activity chair both in school and at home
- Angle adjustable tray is a perfect surface for work, play, communications, nutrition, and sensory activities



Specifications

	Bantam extra small	Bantam small	Bantam medium
Approximate Height Range	28"-40"	36"-54"	4'-5'6"
Maximum User Weight	50 lbs	100 lbs	200 lbs
Seat Depth Range w/Back	7"-12"	11"-16"	15"-20"
Seat to Foot Plate Range	4.5"-16.5"	4.5"-16.5"	7"-18.5"
Seat Height from Floor	19"	19"	22" Planar Seat
Knee Pad Size	2.5"Wx5"H	3.5"Wx5"H	7"Wx9"H
Knee Pad Depth Range	2.5"-8.5"	2.5"-8.5"	3"-8"
Weight of Basic Unit	54 lbs	54 lbs	108 lbs
Frame Footprint	24.5"x36.5"	24.5"x36.5"	24.5"x44.5"
Supine	0-90°	0-90°	0-90°



EasyStand*li*

262 West 1st St • Morton, MN 56270, USA • 800.342.8968 • info@easystand.com • easystand.com

Standing technology should only be used under the guidance of a physician with recommendations for standing program protocol and any medical precautions. Standing programs should be monitored by the attending therapist. AMI maintains a policy of continual product improvement and reserves the right to change features, specifications, and prices without prior notification. Check with AMI for latest information. FORM BANMEDADVAN-072715 REV A Copyright © 2015 Altitude Medical, Inc. All rights reserved. U.S. & international patents pending. Printed in the U.S.A.