



Zen-O™

## Portable Oxygen Concentrator

Give your patients more freedom  
with mobile oxygen on demand



Gas Control Equipment



SIMPLE. PORTABLE. DEPENDABLE.

## MAXIMISE YOUR PATIENTS' INDEPENDENCE AND PEACE OF MIND

You can help many of your oxygen patients to live their lives more independently with the Zen-O™ portable oxygen concentrator from GCE Healthcare.

The mobile, lightweight Zen-O™ will enable most of your patients to live an active lifestyle. Zen-O™ will also help ease any anxieties some might have about being left without supplementary oxygen.



## DESIGNED FOR PATIENTS

We've worked hard to make sure Zen-O™ delivers oxygen to your patients when and where they need it most.

Zen-O™ uses patented technology that quietly adjusts the flow rate of oxygen in response to a patient's breath rate. For example, if your patient climbs some stairs and the device detects the need for more oxygen, it automatically increases the amount delivered to ensure they remain saturated.



## PULSE AND CONTINUOUS FLOW

Zen-O™ gives patients the choice of either pulse or continuous flow rates at up to 2 litres per minute of oxygen produced.

Your business will also benefit, as Zen-O™ enables you to meet the oxygen needs of most of your patients.



SIMPLE. PORTABLE. DEPENDABLE.

## LOW COST TO OXYGEN PROVIDERS

While Zen-O™ delivers freedom, choice and peace of mind to patients, home oxygen providers will find it easy to maintain and so the cost of ownership is low.



## NO WAITING OR RETURNS

When the sieve bed needs to be replaced, there's no need to incur costs by sending it back to GCE Healthcare or an authorised distributor for maintenance. Best of all, your patients don't have to wait for its return.

Zen-O's design means your trained staff can use simple tools to carry out basic maintenance, such as replacing the sieve bed. It just takes a few minutes.

And, after the sieve bed has been successfully replaced, Zen-O™ can be updated using a combination of buttons on the display.

Our global presence and robust supply chain means replacement parts and accessories are readily available.

## FLEXIBILITY TO GO WITH THE FLOW

Weighing just 4.66kg (10.25lbs), the Zen-O™ offers lightweight portability, which means patients can carry on with their lives unrestricted.

With its intuitive buttons and clear LCD display, many of your patients will find the ultra-reliable Zen-O™ simple to operate in all kinds of environments.



## ONE SOLUTION

Zen-O™ comes with everything most patients need to remain mobile, including an unobtrusive carry bag, functional pull cart, handy DC and AC power cords, a user manual and quick reference guide that are easy to understand.

Your device will include a re-chargeable lithium-ion battery for use away from home. Zen-O™ can hold up to two batteries that can last up to eight hours in pulse mode.



SIMPLE. PORTABLE. DEPENDABLE.

## BUILT TO LAST

Manufactured with pride in the UK and to the highest standards, Zen-O™ will perform so quietly, reliably and unobtrusively, your patients will hardly know it's there.

Zen-O™ portable oxygen concentrator is supplied with a three year device warranty or 15,000 hours of total use for extra assurance.

Zen-O™ is built using carefully selected components and accessories. It's designed to withstand everyday use and an internationally recognised testing facility has independently tested and assessed it to extensive and exacting standards.

## MANUFACTURED TO HIGH STANDARDS

Zen-O™ meets the requirements of the;

- European Medical Device Directive
- United States Food and Drug Administration
- United States Federal Aviation Administration



## SPECIFICATIONS

Size (WxDxH):	212 mm × 168 mm × 313 mm (8.3" × 6.6" × 12.3")
Weight:	4.66 kg with one 12 cell battery
Power requirements:	AC adaptor: 100-240V AC(+/- 10%) 50-60 Hz in, 24V DC, 6.25A out DC adaptor: 11.5 - 16V DC in, 19V, 7.9A out
Purity:	87% - 96% at all settings
Maximum oxygen discharge pressure:	20.5 psi
Inspiratory trigger sensitivity:	-0.12cm/H <sub>2</sub> O
Humidity range:	5% to 93% ± 2% non-condensing
Operating altitude:	0 to 9,000 ft relative to sea level (0 km to 2.7 km)
Temperature:	
Operation:	5°C (41°F) and 40°C (104°F)
Storage:	-20°C (-4°F) and 60°C (140°F)
Setting:	Adjustable in 0.5 increments from 1.0 to 6.0 in pulse mode and from 0.5 to 2.0 in continuous mode
Noise level:	38 dB(A) tested according to Prüfmethode 14-1 03/2007 MDS-Hi* 42 dB(A) tested according to ISO 3744*
Alarm types:	Low oxygen purity No breath detected Low battery Service required
Battery duration:	Up to 4 hours with a single battery or 8 hours with 2 batteries at 18 BPM

\*At pulse setting 2.



## EUROPE

### CZECH REPUBLIC

GCE Trade s.r.o.  
Zizkova 381  
CZ-583 01 Chotebor  
Phone: +420 569 661 122  
Fax: +420 569 661 107  
sales.cz@gcegroup.com

### FRANCE

GCE S.A.S.  
70, rue du PUIITS CHARLES  
BP N° 40110  
FR-58403 La Charité-sur-Loire  
Phone: +33/3 86 69 46 00  
Fax: +33/3 86 70 09 15  
sales.fr@gcegroup.com

### GERMANY

GCE GmbH  
Weyherser Weg 8  
36043 Fulda  
Phone: +49 (0)661-8393-0  
Fax: +49 661 8393 25  
sales.de@gcegroup.com

### HUNGARY

GCE Hungária Kft.  
II. Rákóczi Ferenc út 68.  
H-2314 Halásztelek,  
Phone: +36 (24) 521 200  
Fax: +36 (24) 521 201  
sales.hu@gcegroup.com

### ITALY

GCE Mujelli spa  
Via F. Cervi, 11  
37036 San Martino B.A. (VR)  
Phone: +39 045 878 0 525  
Fax: +39 045 878 0 750  
sales.it@gcegroup.com

### POLAND

GCE Sp z o.o.  
ul. Drapińska 12  
03-581 Warszawa  
Phone: +48 22 511 23 57  
Fax: +48 22 677 70 90  
sales.pl@gcegroup.com

### PORTUGAL

GCE Portugal  
Rua do FeiraPark, 50, Piso  
2, sala 8  
PT-4520-632 São João de Vêr  
Phone: +351 256 373 682  
Fax: +351 256 378 260  
sales.pt@gcegroup.com

### ROMANIA

GCE Romania S.R.L  
22, Al.Puskin Street, Ap.1,  
Bucharest 1, 011996  
Phone: +40 21 316 76 72  
Fax: +40 21 316 76 74  
sales.ro@gcegroup.com

### RUSSIA

GCE Krass OOO  
Russian Federation  
129343, Moscow,  
Urzhymkaya Street, no 4  
Phone: +7 495 745 26 99  
Fax: +7 812 323 86 49  
sales.ru@gcegroup.com

### RUSSIA

GCE Krass OOO  
Russian Federation  
194100, Sankt Petersburg  
Kantemirovskaya, 12A  
Phone: +7 812 323 86 21  
Fax: +7 812 323 86 49  
sales.ru@gcegroup.com

### SPAIN

GCE Ibérica, S.L.U.  
Avda. de la Democracia,  
7 - Oficina 311  
ES-28031 Madrid  
Phone: +34 915 711 470  
Fax: +34 915 712 756  
sales.es@gcegroup.com

### SWEDEN

GCE Norden AB  
Box 21044  
Källvattengatan 9  
200 21 Malmö  
Phone: +46-40-388300  
Fax: +46-40-388353  
sales.se@gcegroup.com

## UNITED KINGDOM

GCE Ltd.  
100 Empress Park, Penny Lane,  
Haydock, ST HELENS  
WA11 9DB  
Phone: +44 (0)1942 29 29 50  
Fax: +44 (0)1942 29 29 77  
sales.gb@gcegroup.com

## ASIA

### CHINA

GCE Gas Control Equipment Co., Ltd.  
No.4 Building, 318 Xiao Wan Road  
Fengxian District, Shanghai 201401  
P.R. China  
Phone: +86-21-37198408  
Fax: +86-21-37198617  
sales.cn@gcegroup.com

### INDIA

GCE India Pvt. Ltd.  
59, 1<sup>st</sup> Floor,  
Millers Road, Benson Town,  
Bangalore - 560046, Karnataka  
Phone: +91 802 363 1685  
Fax: +91 802 353 0110  
sales.in@gcegroup.com

### INDIA

GCE India Representation Office  
Offices No. 44 & 45, 2nd Floor,  
Crystal Plaza, Hiranandani  
Complex,  
Sector - 7, Kharghar,  
Navi Mumbai. 410 210,  
Maharashtra, India  
Mobile: +91-9987026546  
sales.in@gcegroup.com

## AMERICAS

### MEXICO

GCE Gas Control Equipment SA de CV  
Miguel de Cervantes Saavedra 193  
Ampliación Granada Miguel Hidalgo  
Distrito Federal 11529, México  
Phone: +52 55 2626 1439  
sales.pa@gcegroup.com

### PANAMA

GCE Latin America  
Po.Box: 0843-01211  
Panamá  
Republica de Panamá  
Phone: +507 317 61 68  
Fax: + 507 317 65 00  
sales.pa@gcegroup.com

### UNITED STATES OF AMERICA & CANADA

Gas Control Equipment Inc.  
1508 Windsor Forest Trail  
Keller, TX 76262  
U.S.A



Gas Control Equipment  
www.gcegroup.com