

# PRODUCT INSERT

## LIFT CHAIR WEIGHT CAPACITIES

Your lift chair was developed for your convenience. See the below specification table for information on weight capacity. Refer to the owner's manual for any additional information concerning your lift chair.

MODEL	WEIGHT CAPACITY
102	250 lbs. (113 kg)
435, 805, 835	300 lbs. (136 kg)
10, 15, 105, 106, 107, 110, 158, 215, SR225	325 lbs. (147 kg)
20, 30, 58, 108, 125S, 125M, 125L, 205, 220, LC225, 250, 310, 325, 358P, 358PW, 358S, 358M, 358L, 421, 450, 450S, 465M, 470S, 470M, 470W, 470LT, 485, 510P, 510S, 510M, 510L, 521, 525S, 525M, 525L, 525PW, 550S, 550M, 550L, 570P, 570PW, 570S, 570M, 570L, 570W, 570T, 571, 575, 576M, 576L, 585, 590, 595, 660, 670, 770S, 770L, 900	375 lbs. (170 kg)
358XL	500 lbs. (227 kg)
358XXL	600 lbs. (272 kg)



**MANDATORY!** Stay within the specified weight capacity of your lift chair. Exceeding the weight capacity voids your lift chair warranty.

**NOTE:** Lift Chair model numbers may contain a prefix (i.e. SR or LC). When using the table above, look for only the numeric portion of your specific model.



Votre fauteuil auto-souleveur a été conçu pour vous faciliter la vie. Voir ci-dessous le tableau des spécifications afin de connaître la limite de poids maximale pour votre fauteuil. Référez-vous au manuel du fauteuil élévateur pour toutes autres informations.

MODELE	CAPACITE POIDS MAXIMUM
102	113 kg (250 lb)
435, 805, 835	136 kg (300 lb)
10, 15, 105, 106, 107, 110, 158, 215, SR225	147 kg (325 lb)
20, 30, 58, 108, 125S, 125M, 125L, 205, 220, LC225, 250, 310, 325, 358P, 358PW, 358S, 358M, 358L, 421, 450, 450S, 465M, 470S, 470M, 470W, 470LT, 485, 510P, 510S, 510M, 510L, 521, 525S, 525M, 525L, 525PW, 550S, 550M, 550L, 570P, 570PW, 570S, 570M, 570L, 570W, 570T, 571, 575, 576M, 576L, 585, 590, 595, 660, 670, 770S, 770L, 900	170 kg (375 lb)
358XL	227 kg (500 lb)
358XXL	272 kg (600 lb)



**OBLIGATOIRE!** N'excédez pas la limite de poids maximale de votre fauteuil auto-souleveur sous peine d'annuler la garantie.

**REMARQUE:** Les numéros de modèle du fauteuil auto-souleveur peuvent contenir un préfixe (exemple: SR ou LC). Lorsque vous utilisez le tableau ci-dessus, regardez uniquement la partie numérique de votre modèle spécifique.

